

Product datasheet for TP322369

GRIK2 (NM_021956) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutamate receptor, ionotropic, kainate 2 (GRIK2), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222369 representing NM_021956 Red=Cloning site Green=Tags(s)

MKIIFPILSNPVFRRTVKLLLCLLWIGYSQGTTHVLRFGGIFEYVESGPMGAEELAFRFAVNTINRNRTL
LPNTTLTYDTQKINLYDSFEASKKACDQLSLGVAFIGPSHSSSANAVQSICNALGVPHIQTRWKHQVSD
NKDSFYVSLYPDFSSLSRAILDVQFFKWKTVT VVYDDSTGLIRLQELIKAPSRYNLRLKIRQLPADTKD
AKPLLKEMKRGKEFHVIFDCSHEMAAGILKQALAMGMMTEYYHYIFTLLDLFALDVEPYRYSVGNMTGFR
ILNTEQVSSIIKWSMERLQAPPKPDSGLLDGFMTTDAALMYDAVHVSVAVQQFPQMTVSSLQCNRH
KPWRFGTRFMSLIKEAHWEGLTGRITFNKTNGLRDTDFLDVISLKEEGLEKIGTWDPASGLNMTESQKKG
PANITDSLSNRSLIVTTILEEPYVLFKKS DKPLYGNDRFEGYCIDLLRELSTILGFTYEIRLVEDGKYGA
QDDANGQWNGMVRELIDHKADLAVAPLAITYVREKVIDFSKPFMTLGISILYRKPNGTNPVGFSLNPLS
PDIWMYVLLAYLGVSCVLFVIARFSPYEWYNPHPCNPDSVWENNFTLLNSFWFGV GALMRQGS ELM PKA
LSTRIVGGIWWFFTLIISSYTANLAAFLTVERMESPID SADDLAKQTKIEYGAVEDGATMTFFKSKIS
TYDKMWAFMSSRRQSVLVKSNEEGIQRVLTSDYAFLMESTTIEFVTQRNCNLTQIGGLIDSKGYGVGTPM
GSPYRDKITIAILQLQEEGKLHMMKEKWWRGNGCPEEESKEASALGVQNI GGIFIVLAAGLVLSV FVAVG
EFLYKSKKNAQLEKRSFCSAMVEELRMSLKCQRRLKHKPQAPVIVKTEEVINMHTFNDRRLPGKETMA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	98.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_068775
Locus ID:	2898
UniProt ID:	Q13002 , Q8IY40 , A8K0H7
RefSeq Size:	3322
Cytogenetics:	6q16.3
RefSeq ORF:	2724
Synonyms:	EAA4; GLR6; GluK2; GLUK6; GLUR6; MRT6
Summary:	Glutamate receptors are the predominant excitatory neurotransmitter receptors in the mammalian brain and are activated in a variety of normal neurophysiologic processes. This gene product belongs to the kainate family of glutamate receptors, which are composed of four subunits and function as ligand-activated ion channels. The subunit encoded by this gene is subject to RNA editing at multiple sites within the first and second transmembrane domains, which is thought to alter the structure and function of the receptor complex. Alternatively spliced transcript variants encoding different isoforms have also been described for this gene. Mutations in this gene have been associated with autosomal recessive cognitive disability. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Ion Channels: Glutamate Receptors, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction

Product images:



Coomassie blue staining of purified GRIK2 protein (Cat# TP322369). The protein was produced from HEK293T cells transfected with GRIK2 cDNA clone (Cat# [RC222369]) using MegaTran 2.0 (Cat# [TT210002]).