

Product datasheet for TP322395

DGKA (NM_001345) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens diacylglycerol kinase, alpha 80kDa (DGKA), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222395 representing NM_001345 Red =Cloning site Green =Tags(s)

MAKERGLISPSDFAQLQKYMESTKKVSDVLKLFEDGEMAKYVQGDAIGYEGFQQFLKIYLEVDNVPRL
SLALFQSFETGHCLNETNVTKDWCCLNDVSCYFSLLEGGRPEDKLEFTFKLYDTRNGILDSSEVDKIIL
QMMRVAEYLDWDVSELRPILQEMMKEIDYDGSQSVSQAEWVRAGATTVPLLVLGLEM TLKDDGQHM
WRP
KRFPRPVYCNLCESIGLGKQGLSCNLCKYTVHDQCAMKALPCEVSTYAKSRKDIGVQSHVWVRGGCESG
RCDRCQKKIRIYHSLTGLHCVWCHLEIHDDCLQAVGHECDCGLLRDHILPPSSIYPSVLASGPDRKNSKT
SQKTMDLNLSTSEALRIDPVPNTHPLLVFVNPKSGGKQGQVRLWKFQYILNPRQVFNLLKDGPEIGLRL
FKDVPDSRILVCGGDGTGWILETIDKANLPVLPVAVLPLGTGNDLARCLRWGGGYEQNLAKILKDLE
MSKVVHMDRWSVEVIPQQTEEKSDPVPFQIINNYFSIGVDASIAHRFHIMREKYPEKFNSRMKNKLWYFE
FATSESIFFSTCKKLEESLTVEICGKPLDLSNLSLEGIAVLNIPSMHGGSNLWGDTRRPHGDIYGINQALG
ATAKVITDPDILKTCVPLSDKRLEVVGLEGAIEMGQIYTKLKNAGRRLAKCSEITFHTTKLPMQIDGE
PWMQTPCTIKITHKNQMPMLMGPPRSTNFFGFLS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	82.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001336
Locus ID:	1606
UniProt ID:	P23743
RefSeq Size:	2756
Cytogenetics:	12q13.2
RefSeq ORF:	2205
Synonyms:	DAGK; DAGK1; DGK-alpha
Summary:	The protein encoded by this gene belongs to the eukaryotic diacylglycerol kinase family. It acts as a modulator that competes with protein kinase C for the second messenger diacylglycerol in intracellular signaling pathways. It also plays an important role in the resynthesis of phosphatidylinositols and phosphorylating diacylglycerol to phosphatidic acid. Several transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Apr 2017]
Protein Families:	Druggable Genome
Protein Pathways:	Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways, Phosphatidylinositol signaling system

Product images:



Coomassie blue staining of purified DGKA protein (Cat# TP322395). The protein was produced from HEK293T cells transfected with DGKA cDNA clone (Cat# [RC222395]) using MegaTran 2.0 (Cat# [TT210002]).