

Product datasheet for TP322395M

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

DGKA (NM_001345) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Homo sapiens diacylglycerol kinase, alpha 80kDa (DGKA),

transcript variant 3, 100 µg

Species: Human Expression Host: HEK293T

Expression cDNA >RC222395 representing NM_001345
Clone or AA Red=Cloning site Green=Tags(s)

Sequence:

MAKERGLISPSDFAQLQKYMEYSTKKVSDVLKLFEDGEMAKYVQGDAIGYEGFQQFLKIYLEVDNVPRHL SLALFQSFETGHCLNETNVTKDVVCLNDVSCYFSLLEGGRPEDKLEFTFKLYDTDRNGILDSSEVDKIIL QMMRVAEYLDWDVSELRPILQEMMKEIDYDGSGSVSQAEWVRAGATTVPLLVLLGLEMTLKDDGQHMWRP

KRFPRPVYCNLCESSIGLGKQGLSCNLCKYTVHDQCAMKALPCEVSTYAKSRKDIGVQSHVWVRGGCESG RCDRCQKKIRIYHSLTGLHCVWCHLEIHDDCLQAVGHECDCGLLRDHILPPSSIYPSVLASGPDRKNSKT SQKTMDDLNLSTSEALRIDPVPNTHPLLVFVNPKSGGKQGQRVLWKFQYILNPRQVFNLLKDGPEIGLRL FKDVPDSRILVCGGDGTVGWILETIDKANLPVLPPVAVLPLGTGNDLARCLRWGGGYEGQNLAKILKDLE MSKVVHMDRWSVEVIPQQTEEKSDPVPFQIINNYFSIGVDASIAHRFHIMREKYPEKFNSRMKNKLWYFE FATSESIFSTCKKLEESLTVEICGKPLDLSNLSLEGIAVLNIPSMHGGSNLWGDTRRPHGDIYGINQALG ATAKVITDPDILKTCVPDLSDKRLEVVGLEGAIEMGQIYTKLKNAGRRLAKCSEITFHTTKTLPMQIDGE

 ${\tt PWMQTPCTIKITHKNQMPMLMGPPPRSTNFFGFLS}$

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 82.5 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.



DGKA (NM_001345) Human Recombinant Protein - TP322395M

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001336

Locus ID: 1606

UniProt ID: <u>P23743</u>, <u>A0A024RB23</u>

RefSeq Size: 2756
Cytogenetics: 12q13.2
RefSeq ORF: 2205

Synonyms: DAGK; DAGK1; DGK-alpha

Summary: The protein encoded by this gene belongs to the eukaryotic diacylglycerol kinase family. It acts as

a modulator that competes with protein kinase C for the second messenger diacylglycerol in

intracellular signaling pathways. It also plays an important role in the resynthesis of

phosphatidylinositols and phosphorylating diacylglycerol to phosphatidic acid. Several transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Apr


2017]

Protein Families: Druggable Genome

Protein Pathways: Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways,

Phosphatidylinositol signaling system

Product images:



Coomassie blue staining of purified DGKA protein (Cat# [TP322395]). The protein was produced from HEK293T cells transfected with DGKA cDNA clone (Cat# [RC222395]) using MegaTran 2.0 (Cat# [TT210002]).