

Product datasheet for TP322532L

SLC26A4 (NM_000441) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens solute carrier family 26, member 4 (SLC26A4), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222532 representing NM_000441 Red =Cloning site Green =Tags(s)

MAAPGGRSEPPQLPEYSCSYMVS RPVYSELAFQQQHERRLQERKTLRESLAKCCSCSRKRAFGVLKTLVP
 ILEWLPKYRVKEWLLSDVISGVSTGLVATLQGMAYALLA AVPVGYGLYSAFFPILTYFIFGTSRHISVGP
 FPVSLMVGSVVLSMAPDEHFLVSSSNGTVLNTTMDTAARDTARVLIASALTLLVGIIQLIFGG LQIGF
 IVRYLADPLVGGFTTAAAFQVLVSQ LKIVLNVSTKKNYNGVLSIIYTLVEIFQNIGDNTLADFTAGLLTIV
 VCM AVKELNDRFRHKIPVPIPIEVIVTIIATAISYGANLEKNYNAGIVKSIPRGFLPPELPPVSLFSEML
 AASFSAVVAYAI AVSVGKVYATKYDYTIDGNQEFIAFGISNIFSGFFSCFVATTALSRTAVQESTGGKT
 QVAGIISAAIVMIAIALGKLEPLQKSVLAAVVIANLKGMMQLCDIPRLWRQNKIDAVI WVFTCIVSI
 ILGLDLGLLAGLIFGLLT VVLRVQFPSWNGLSIPSTDIYKSTKKNYKNIIEPQGVKILRFSSPIFYGNVD
 GFKKCIKSTVGFD AIRVYNKRLKALRKIQKLIKSGQLRATKNGIISDAVSTNNAFEPDEDIEDLEELDIP
 TKEIEIQVDWNS ELPVKVNVPKVPIHSLVLD CGAISFLDVVGVRSLRVIVKEFQRIDVNVYFASLQDYVI
 EKLEQCGFFDDNIRKDTFFLTVHDAILYLQNQVKSQEGQGSILETITLIQDCKDTLELIETELTEEELDV
 QDEAMRTLAS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	85.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000432
Locus ID:	5172
UniProt ID:	O43511
RefSeq Size:	4930
Cytogenetics:	7q22.3
RefSeq ORF:	2340
Synonyms:	DFNB4; EVA; PDS; TDH2B
Summary:	Mutations in this gene are associated with Pendred syndrome, the most common form of syndromic deafness, an autosomal-recessive disease. It is highly homologous to the SLC26A3 gene; they have similar genomic structures and this gene is located 3' of the SLC26A3 gene. The encoded protein has homology to sulfate transporters. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified SLC26A4 protein (Cat# [TP322532]). The protein was produced from HEK293T cells transfected with SLC26A4 cDNA clone (Cat# [RC222532]) using MegaTran 2.0 (Cat# [TT210002]).