

Product datasheet for TP322614

CD41 (ITGA2B) (NM_000419) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human integrin, alpha 2b (platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41) (ITGA2B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222614 protein sequence Red =Cloning site Green =Tags(s)

MARALCPLQALWLLEWLLLLGPCAAPPAAWALNLDPVQLTFYAGPNGSQFGFSLDFHKDSHGRVAIVWG
A
PRTLGPSQEETGGVFLCPWRAEGGQCPSLLFLDRDETRNVGSQTLQTFKARQGLGASVSWSDVIVACAP
WQHWNVLEKTEEAETKTPVGSFLAQPEGRRAEYSPCRGNTLSRIYVENDFSWDKRYCEAGFSSVVTQAG
ELVLGAPGGYYFLGLLAQAPVADIFSSYRPGILLWHVSSQSLSFDSSNPEYFDGYWGYSVAVGEFDGDLN
TTEYVVGAPTWSWTLGAVEILDSYYQRLHRLRGEQMASYFGHSVAVTDVNGDGRHDLLVGAPLYMESRA
D
RKLAEVGRVYLFQPRGPHALGAPSLLLTGTQLYGRFGSAIAPLGDLDLDRDGYNDIAVAAPYGGPSGRGQV
LVFLGQSEGLRSRPSQVLDSPFPTGSFGLRGAVIDDNGYPDLIVGAYGANQVAVYRAQPVKASVQ
LLVQDSLNPVAVKSCVLPQTKTPVSCFNIQMCVGTGHNIPQKLSLNAELQLDRQKPRQGRRVLLLSGQQ
A
GTTLNLDLGGKHSPICHTTMAFLRDEADFRDKLSPIVLSLNVSLPPEAGMAPAVWLHGDTHTVQEQRTRIV
LDCGEDDVCVPQLQLTASVTGSPLLVGADNVLELQMDAANEGEGAYEAEAVHLPQGAHYMRALSNEVEG
F
ERLICNQKKENETRVLCELGNPMKKNAQIGIAMLVSVGNLEEAGESVSFQLQIRSKNSQNPNSKIVLLD
VPVRAEAQVELRGNSFPASLVAAEEGEREQNSLDSWGPKEHTYELHNNPGTVNGLHLSIHLPGQSQ
P
SDLLYILDIQPQGGQLQCFPPVNPVKVDWGLPIPSPIHPAHHKRDRRQIFLPEPEQPSRLQDPVLVS
CDSAPCTVVQCDLQEMARGQRAMVTVLAFWLPLSLYQRPLDQFVLQSHAWFNVSSLPYAVPPLSLPRGE
A
QVWTQLLRALERAIPWVWLVGVLGLLLLLILVLMWVKVGFKRNRPPEEDDEEGE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK



[View online »](#)

Predicted MW:	110 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000410
Locus ID:	3674
UniProt ID:	P08514
RefSeq Size:	3334
Cytogenetics:	17q21.31
RefSeq ORF:	3117
Synonyms:	BDPLT2; BDPLT16; CD41; CD41B; GP2B; GPIIb; GT; GT1; GTA; HPA3; PPP1R93
Summary:	This gene encodes a member of the integrin alpha chain family of proteins. The encoded preproprotein is proteolytically processed to generate light and heavy chains that associate through disulfide linkages to form a subunit of the alpha-IIb/beta-3 integrin cell adhesion receptor. This receptor plays a crucial role in the blood coagulation system, by mediating platelet aggregation. Mutations in this gene are associated with platelet-type bleeding disorders, which are characterized by a failure of platelet aggregation, including Glanzmann thrombasthenia. [provided by RefSeq, Jan 2016]
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hematopoietic cell lineage, Hypertrophic cardiomyopathy (HCM), Pathways in cancer, Regulation of actin cytoskeleton, Small cell lung cancer

Product images:

Coomassie blue staining of purified ITGA2B protein (Cat# TP322614). The protein was produced from HEK293T cells transfected with ITGA2B cDNA clone (Cat# [RC222614]) using MegaTran 2.0 (Cat# [TT210002]).