

Product datasheet for TP322687L

ZNF148 (NM_021964) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens zinc finger protein 148 (ZNF148), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222687 representing NM_021964 Red=Cloning site Green=Tags(s)

MNIDDKLEGLFLKCGGIDEMQSSRTMVVMGGVSGQSTVSGELQDSVLQDRSMPHQEILAADEVLQESEMR
QQDMISHDELMVHEETVKNDEEQMETHERLPQGLQYALNPISVKQEITFTDVSEQLMRDQKQIREPVDL
QKKKKRKRSPAKILTINEDGSLGLKTPKSHVCEHCNAAFRTNYHLQRHVFIHTGEKPFQCSQCDMRFIQ
KYLQRHEKIHTGEKPFRCDECGMRFIQKYHMERHKRTHSGEKPYQCEYCLQYFSRTDRVLKHKRMCHEN
HDKKLNRCIAKGGLLTSEEDSGFSTSPKDNSLPKRRQKTEKSSGMDKESALDKSDLKKDKNDYLPPLYS
SSTKVKDEYMVAEYAVEMPHSSVGGSHLEDASGEIHPPKLVLLKINSKRSLKQPLEQNQTISPLSTYEES
KVSXYAFELVDKQALLDSEGNADIDQVDNLQEGPSKPVHSSSTNYDDAMQFLKKRYLQAASNSREYALN
VGTIASQPSVTQAAVASVIDESTTASILESQALNVEIKSNHDKNVIPDEVLTLLDHYSHKANGQHEISF
SVADTEVTSSISINSSEVPEVTPSENVGSSSQASSSDKANMLQEYSKFLQALDRTSQNDAYLNSPSLNF
VTDNQTLPNQPAFSSIDKQVYATMPINSFRSGMNSPLRTPDKSHFGLIVGDSQHSFPFSGDETNHASAT
STQDFLDQVTSQKAEAPVHQAYQMSSFEQPFPRAPYHGSRAGIATQFSTANGQVNLRGPPTSAEFSEFP
LVNVNDNRAGMTSSPDATTGQTFG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	88.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_068799
Locus ID:	7707
UniProt ID:	Q9UQR1
RefSeq Size:	3032
Cytogenetics:	3q21.2
RefSeq ORF:	2382
Synonyms:	BERF-1; BFCOL1; GDACCF; HT-BETA; pHZ-52; ZBP-89; ZFP148
Summary:	The protein encoded by this gene is a member of the Kruppel family of zinc finger DNA binding proteins. The encoded protein activates transcription of the T-cell receptor and intestinal alkaline phosphatase genes but represses transcription of the ornithine decarboxylase, vimentin, gastrin, stomelysin, and enolase genes. Increased expression of this gene results in decreased patient survival rates from colorectal cancer, while mutations in this gene have been associated with global developmental delay, hypoplastic corpus callosum, and dysmorphic facies. [provided by RefSeq, Feb 2017]
Protein Families:	Transcription Factors

Product images:



Coomassie blue staining of purified ZNF148 protein (Cat# [TP322687]). The protein was produced from HEK293T cells transfected with ZNF148 cDNA clone (Cat# [RC222687]) using MegaTran 2.0 (Cat# [TT210002]).