

Product datasheet for **TP322714**

THYN1 (NM_199297) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human thymocyte nuclear protein 1 (THYN1), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222714 protein sequence Red =Cloning site Green =Tags(s)
	 MSRPRKRLAGTSGSDKGLSGKRTKTENSGEALAKVEDSNPQKTSATKNCLKNLSSHWMKSEPEERLEKG VDVKFSIEDLKAQPKQTTTCWDGVRNYQARNFLRAMKLGEEAFFYHSNCKEPIAGLMKIVKEAYPDHTQF EKNNPHYDPSSKEDNPKWSMKSLILF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	18.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_954994
Locus ID:	29087
UniProt ID:	Q9P016



[View online »](#)

RefSeq Size: 862


Cytogenetics: 11q25

RefSeq ORF: 498

Synonyms: HSPC144; MDS012; MY105; THY28; THY28KD

Summary: This gene encodes a protein that is highly conserved among vertebrates and plant species and may be involved in the induction of apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified THYN1 protein (Cat# TP322714). The protein was produced from HEK293T cells transfected with THYN1 cDNA clone (Cat# [RC222714]) using MegaTran 2.0 (Cat# [TT210002]).