

## Product datasheet for **TP322724L**

### **CELF2 (NM\_001083591) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Homo sapiens CUG triplet repeat, RNA binding protein 2 (CUGBP2), transcript variant 4, 1 mg
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>RC222724 representing NM_001083591 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MNGALDHSQDPDPAIKMFVGGQIPRSWSEKELKELFEPYGA VYQINVLDRSQNPPQSKGCCFVTFYTRK AALEAQNALHNIKTLPGMHHP IQMKPADSEKSNAVEDRKLFIGMVSKKCNENDIRVMFSPFGQIEECRIL RGPDGLSRGCAFVTFSTRAMAQNAIKAMHQSQTMEGCSSPIVVKFADTQKDKERRLQQQLAQMQQLNT ATWGNLTGLGGLTPQYLALLQQATSSSNLGA FSGIQQMAGMNALQLQNLATLAAAAAAAAQTSATSTNANP LSTSSALGALTSPVAASTPNSTAGAAMNSL TSLGTLQGLAGATVGLNNINALAGTINSMAALNGGLGAT GLTNGTAGTMDALTQAYSIGIQYAAAAALPTLYS QLLQQSAAGSQKEGPEGANLFIYHLPQEFGDQDIL QMFMPPFGNVISAKVFIDKQTNLSKCFGFVSYD NPVSAQAAIQAMNGFQIGMKRLKVQLKRSKNDSKPY  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
<b>Tag:</b>	C-Myc/DDK
<b>Predicted MW:</b>	51.9 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_001077060](#)

Locus ID: 10659

UniProt ID: [Q95319](#)

RefSeq Size: 8053


Cytogenetics: 10p14

RefSeq ORF: 1464

Synonyms: BRUNOL3; CELF-2; CUG-BP2; CUGBP2; ETR-3; ETR3; NAPOR

**Summary:** Members of the CELF/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing, and translation. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified CELF2 protein (Cat# [TP322724]). The protein was produced from HEK293T cells transfected with CELF2 cDNA clone (Cat# [RC222724]) using MegaTran 2.0 (Cat# [TT210002]).