

Product datasheet for **TP322797M**

SGK196 (POMK) (NM_032237) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human protein kinase-like protein SgK196 (SGK196), 100 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC222797 protein sequence
Red=Cloning site **Green**=Tags(s)

MEKQPQNSRRGLAPREVPPAVGLLLIMALMNTLLYLCLDHFFIAPRQSTVDPTHCPYGHFRIGQMKNCSPLWLSCEELRTEVRQLKRVGEGAVKRVFLSEWKEHKVALSQLTSLEMKDDFLHGLQMLKSLQGTHVTTLLGYCEDDNTMLTEYHPLGSLNLEETLNLSKYQNVNTWQHRLELAMDYVSIINYLLHHSVPVTRVMCDNSNDLPKTLNQYLLTSNFSILANDLDALPLVNHSSGMLVKCGHRELHGDFVAPEQLWPYGEDVPPFHDDLMPYDEKIDIWKIPDISSFLLGHIEGSDMVRFHLDIHKACKSQTTPSERPTAQDVLETYQKVLDTLRDAMMSQAREML

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 39.9 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_115613](#)

Locus ID: 84197



[View online »](#)

UniProt ID: [Q9H5K3](#)
RefSeq Size: 1623
Cytogenetics: 8p11.21
RefSeq ORF: 1050
Synonyms: MDDGA12; MDDGC12; SGK196

Summary: This gene encodes a protein that may be involved in the presentation of the laminin-binding O-linked carbohydrate chain of alpha-dystroglycan (a-DG), which forms transmembrane linkages between the extracellular matrix and the exoskeleton. Some pathogens use this O-linked carbohydrate unit for host entry. Loss of function compound heterozygous mutations in this gene were found in a human patient affected by the Walker-Warburg syndrome (WWS) phenotype. Mice lacking this gene contain misplaced neurons (heterotopia) in some regions of the brain, possibly from defects in neuronal migration. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, May 2013]

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Product images:



Coomassie blue staining of purified POMK protein (Cat# [TP322797]). The protein was produced from HEK293T cells transfected with POMK cDNA clone (Cat# [RC222797]) using MegaTran 2.0 (Cat# [TT210002]).