

Product datasheet for **TP322820L**

TLR4 (NM_138554) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human toll-like receptor 4 (TLR4), transcript variant 1, 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC222820 representing NM_138554
Red=Cloning site **Green**=Tags(s)

MSASRLAGTLIPAMAFSLCVRPESWEPCEVWPNIYQCMELNFYKIPDNLPFSTKNLDLSFNPLRHLGS
YSFFSPPELQVLDLSRCEIQTIEDGAYQSLSHLSTLILTGNPIQSLALGAFSGLSSLQKLVAVETNLASL
ENFPIGHLKTLKELNVAHNLIQSFKLPEYFSNLTNLEHLDLSSNKIQSIYCTDLRVLHQMPLLNLSDLS
LNPMNFIQPGAFKEIRLHKLTLRNNFDSLNVMTKCIQGLAGLEVHRLVLGEFRNEGNLEKFDKSALEGLC
NLTIIEFRLAYLDYLLDDIIDLFNCLTNVSSFSLSVSTIERVKDFSYNFGWQHLELVNCKFGQPPTLKLK
SLKRLTFTSNKGGNAFSEVDLPSLEFLDLSRNGLSFKGCCSQSDFGTTSLKYLDLSFNGVITMSSNFLGL
EQLEHLDFQHSNLKQMSEFSVFLSLRNLIYLDISHTHTRVAFNGIFNGLSSLEVLKMAGNSFQENFLPDI
FTELRNLTFLDLSQCQLEQLSPTAFNSLSSLQVLNMSHNNFFSLDTPYKCLNSLQVLDYSLNHIMTSKK
QELQHFPSSLAFLNLTQNDFACTCEHQSFQWIKDQRQLLVEVERMECATPSDKQGMPVLSLNITCQMNK
TIIGVSVLSVLVSVWAVLVYKFYFHLMLLAGCIKYGRGENIYDAFVIYSSQDEDWVRNELVKNLEEGVP
PFQLCLHYRDFIPGVAIAANIIHEGFHKSARKVIVVVSQHFIQSRWCIFEYEAQTWQFLSSRAGIIFVL
QKVEKTLRQQVELYRLLSRNTYLEWEDSVLGRHIFWRRLRKALLDGKSWNPEGTVGTGCNWQEATS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag: C-Myc/DDK
Predicted MW: 93.2 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_612564
Locus ID:	7099
UniProt ID:	O00206
RefSeq Size:	5503
Cytogenetics:	9q33.1
RefSeq ORF:	2514
Synonyms:	ARMD10; CD284; TLR-4; TOLL
Summary:	<p>The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. In silico studies have found a particularly strong binding of surface TLR4 with the spike protein of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), the causative agent of Coronavirus disease-2019 (COVID-19). This receptor has also been implicated in signal transduction events induced by lipopolysaccharide (LPS) found in most gram-negative bacteria. Mutations in this gene have been associated with differences in LPS responsiveness, and with susceptibility to age-related macular degeneration. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2020]</p>
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Pathogenic Escherichia coli infection, Toll-like receptor signaling pathway

Product images:



Coomassie blue staining of purified TLR4 protein (Cat# [TP322820]). The protein was produced from HEK293T cells transfected with TLR4 cDNA clone (Cat# [RC222820]) using MegaTran 2.0 (Cat# [TT210002]).