

Product datasheet for **TP322863**

C5orf49 (NM_001089584) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Homo sapiens chromosome 5 open reading frame 49 (C5orf49), 20 µg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC222863 representing NM_001089584 Red =Cloning site Green =Tags(s) |
| | MEDDEEETTASTLRGKPRPPVSAQSAFSYIPRRRLDPKEHSYYRPARTGIISLYDCIFKRRLDYDQKL HRDDREHAKSLGLHVNEEQERPVGVLTSVYGKRINQPIEPLNRDFGRANHVQADFYRKNDIPSLKEPG FGHIAPS TRTRPLEQKLISEEDLAANDILDYKDDDDKV |
| Tag: | C-Myc/DDK |
| Predicted MW: | 16.8 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_001083053 |
| Locus ID: | 134121 |
| UniProt ID: | A4QMS7 |




[View online »](#)

RefSeq Size: 1026

Cytogenetics: 5p15.31

RefSeq ORF: 441

Product images:



Coomassie blue staining of purified C5orf49 protein (Cat# TP322863). The protein was produced from HEK293T cells transfected with C5orf49 cDNA clone (Cat# [RC222863]) using MegaTran 2.0 (Cat# [TT210002]).