

Product datasheet for **TP322903M**

HOPX (NM_139212) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens HOP homeobox (HOPX), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222903 representing NM_139212 Red =Cloning site Green =Tags(s)
	MSAETASGPTEDQVEILEYNFNKVDKHPDSTTLCLIAAEAGLSEEETQKWFQRLAKWRRSEGLPSECRS VTD
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	8.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_631958
Locus ID:	84525
UniProt ID:	Q9BPY8
RefSeq Size:	1116



[View online »](#)

Cytogenetics: 4q12


RefSeq ORF: 219

Synonyms: CAMEO; HOD; HOP; LAGY; NECC1; OB1; SMAP31; TOTO

Summary: The protein encoded by this gene is a homeodomain protein that lacks certain conserved residues required for DNA binding. It was reported that choriocarcinoma cell lines and tissues failed to express this gene, which suggested the possible involvement of this gene in malignant conversion of placental trophoblasts. Studies in mice suggest that this protein may interact with serum response factor (SRF) and modulate SRF-dependent cardiac-specific gene expression and cardiac development. Multiple alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Feb 2009]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified HOPX protein (Cat# [TP322903]). The protein was produced from HEK293T cells transfected with HOPX cDNA clone (Cat# [RC222903]) using MegaTran 2.0 (Cat# [TT210002]).