

Product datasheet for **TP323130M**

CCDC9B (NM_207380) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 15 open reading frame 52 (C15orf52), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223130 representing NM_207380 Red =Cloning site Green =Tags(s) MISCAEQRSRQGEAGRGPAVAPAFPLWLPRGCSGILSVPAVAMHSAGTPRAESPMSRQEKDAELDRRI VALRKKNQALLRRYQEIQEDRRQAEQGGMVTPALLQPDGLTVTISQVPGEKRVVSRNWARGTCGPRV T NEMLEDEDAEDHGGTFCLGELVELAVTMENKAEGKRIVSEKPTRARNQGIEGSPGGRVTRSPPTQVAISS DSARKGSWEPWSRPVGEPPEAGWDYAQWKQEREQIDLARLARHRDAQGDWRRPWDLDKAKSTLQDCS QLR GEGPARAGSRRGPRSHQKLQPPPLPDGKGRGGQASRPSVAPATGSKARGKERLTGRARRWDMKEDKEE L EGQEGSQSTRETPSEEEQAQKQSGMEQGRLGSAPAASPALASPEGPKGESVASTASSVPCSPQEPDLAPL DLSLGGAGIPGRESGCVLGLRPGAQESPVSWPEGSKQQPLGWSNHQAELEVQTCPEPQRGAGLPEPGE D RSGKSGAQQLAPRSRPTRGGSQRSRGTAGVRRRTGRPGPAGRC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	57.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_997263</u>
Locus ID:	388115
UniProt ID:	<u>Q6ZUT6</u>
RefSeq Size:	5344
Cytogenetics:	15q15.1
RefSeq ORF:	1602
Synonyms:	C15orf52

Product images:



Coomassie blue staining of purified CCDC9B protein (Cat# [TP323130]). The protein was produced from HEK293T cells transfected with CCDC9B cDNA clone (Cat# [RC223130]) using MegaTran 2.0 (Cat# [TT210002]).