

## Product datasheet for TP323206L

### DNMT3B (NM\_006892) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human DNA (cytosine-5-)-methyltransferase 3 beta (DNMT3B), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223206 representing NM_006892 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MKGDTRHLNGEEDAGGREDSILVNGACSDQSSDSPPILEAIRTPAIRGRRSSSRLSKREVSSLLSYTQDL  
TGDGDGEDGDGSDTPVMPKLFRETRTRSESPAVRTRNNNSVSSRERHRPSRSTRGRQGRNHVDESPVEF  
PATRSLRRRATASAGTPWPSPPSSYLTDLTDDTEDTHGTPQSSSTPYARLAQDSQQGGMESPVQVADSG  
DGDSEYQDQKEFGIGDLVWVGKIKGFSWVWPAMVVSWKATSKRQAMSGMRWVQWFGDGKFSEVSADKLVAL  
GLFSQHFNLATFNKLVSYRKAMYHALEKARVRAGKTFPSSPGDSLEDQLKPMLEWAHGGFKPTGIEGLKP  
NNTQPVVNKSQVRRRAGSRKLESRKYENKTRRRRTADDSATSDYCPAPKRLKTCYNNNGKDRGDEDQSREQM  
ASDVANNKSSLEDGCLSCGRKNPVSFHPLFEGGLCQTCRDRFLELFYMYDDDGYSYCTVCCEGRELLLC  
SNTSCRCFCVCELEVLVGTGTAAEAKLQEPWSCYMCLPQRCHGVLRRRKDWNVRLQAFFTSDTGLEEYA  
PKLYPAIPAARRRPIRVLSLFDGIATGYLVLVELGKIKVGVASEVCEESIAVGTVKHEGNIKYVNDVRN  
ITKKNIEEWGPFDLVIGGSPCNDLSNVNPARKGLYEGTGRLFFEFYHLLNYSRPEKGDPRPFWMFENVV  
AMKVGDKRDISRFLECNPVMIDAIVSAAHRARYFWGNLPGMNRPVIAKNDKLELQDCLEYNRIAKLKK  
VQTITTKSNSIKQGNQLFPVVMNGKEDVLWCTELERIFGFPVHYTDVSNMGRGARQKLLGRSWSVPVIR  
HLFAPLKDYFACE

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	95.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	20 mM PB, pH 7.6, 10 mM NaCl
Bioactivity:	Surface Plasmon Resonance (SPR) (PMID: <a href="#">27468168</a> )



[View online »](#)

<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_008823</a>
<b>Locus ID:</b>	1789
<b>UniProt ID:</b>	<a href="#">Q9UBC3</a>
<b>RefSeq Size:</b>	4353
<b>Cytogenetics:</b>	20q11.21
<b>RefSeq ORF:</b>	2559
<b>Synonyms:</b>	ICF; ICF1; M.Hsa11B
<b>Summary:</b>	CpG methylation is an epigenetic modification that is important for embryonic development, imprinting, and X-chromosome inactivation. Studies in mice have demonstrated that DNA methylation is required for mammalian development. This gene encodes a DNA methyltransferase which is thought to function in de novo methylation, rather than maintenance methylation. The protein localizes primarily to the nucleus and its expression is developmentally regulated. Mutations in this gene cause the immunodeficiency-centromeric instability-facial anomalies (ICF) syndrome. Eight alternatively spliced transcript variants have been described. The full length sequences of variants 4 and 5 have not been determined. [provided by RefSeq, May 2011]
<b>Protein Families:</b>	Druggable Genome, Embryonic stem cells, Induced pluripotent stem cells, Stem cell - Pluripotency
<b>Protein Pathways:</b>	Cysteine and methionine metabolism, Metabolic pathways

### Product images:



Coomassie blue staining of purified DNMT3B protein (Cat# [TP323206]). The protein was produced from HEK293T cells transfected with DNMT3B cDNA clone (Cat# [RC223206]) using MegaTran 2.0 (Cat# [TT210002]).

