

Product datasheet for TP323253

ANKDD1A (NM_182703) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ankyrin repeat and death domain containing 1A (ANKDD1A), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC223253 representing NM_182703
Clone or AA	Red=Cloning site Green=Tags(s)
Sequence:	MQEELAWETDGLLPLERQLHEAARQNNVGRMQELIGRRVNTRARNHVGRVALHWAAGAGHEQAVRLLLEH EAAVDEEDAVGALTEARLCFGMNALLLSAWFGHLRLILQILVNSGAKIHCESKDGLTLLHCAAQKGHPVL AFIMEDLEDVALDHVDKLGRTAFHRAAEHGQLDALDFLVGSGCDHNVKDKEGNTALHLAAGRGHMAVLQR LVDIGLDLEEQNAEGLTALHSAAGGSHPDQVQLLLRAGSTVNALTQKNLSCLHYAALSGSEDSRVLIHA GGCANVVDHQGASPLHLAVRHNFPALVRLLINSDDVNAVDNRQQTPLHLAAEHAWQDIADMILLIAGVDL NLRDKQGKTALAVAVRSNHVSLVDMIIKADRFYRWEKDHPSPSGKSLSFQDHRQETQQLRSVLWRLAS RYLQPREWKKLAYSWEFTEAHVDAIEQQWTGTRSYQEHGHRMLLIWLHGVATAGENPSKALFEGLVAIGR RDLAGWSTMARSQLTATSASRVQMILVPQPPE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	57.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_874362</u>
Locus ID:	348094
UniProt ID:	<u>Q495B1</u>
RefSeq Size:	3126
Cytogenetics:	15q22.31
RefSeq ORF:	1566

Product images:


Coomassie blue staining of purified ANKDD1A protein (Cat# TP323253). The protein was produced from HEK293T cells transfected with ANKDD1A cDNA clone (Cat# [RC223253]) using MegaTran 2.0 (Cat# [TT210002]).