

## Product datasheet for **TP323273L**

### **NKIRAS2 (NM\_001001349) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NFKB inhibitor interacting Ras-like 2 (NKIRAS2), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223273 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MGKSCKWVCGQASVVGKTSILEQLLYGNHVVGSEMIETQEDIYVGSJETDRGVREQVRFYDTRGLRDGAE LPRHCFSCTDGYVLVYSTDSRESFQRVELLKKEIDKSKDKKEVTIVVLGNKCDLQEQRVDPDVAQHWAK SEKVKLWEVSVADRRSLLEPFVYLASKMTQPQSKSAFPLSRKNKSGSLDG  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	21.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001001349</a></u>
Locus ID:	28511
UniProt ID:	<u><a href="#">Q9NYR9</a></u> , <u><a href="#">A0A024R1Z4</a></u>



[View online »](#)

RefSeq Size: 2389


Cytogenetics: 17q21.2

RefSeq ORF: 573

Synonyms: kappaB-Ras2; KBRAS2

**Summary:** Atypical Ras-like protein that acts as a potent regulator of NF-kappa-B activity by preventing the degradation of NF-kappa-B inhibitor beta (NFKBIB) by most signals, explaining why NFKBIB is more resistant to degradation. May act by blocking phosphorylation of NFKBIB and nuclear localization of p65/RELA NF-kappa-B subunit. It is unclear whether it acts as a GTPase. Both GTP- and GDP-bound forms block phosphorylation of NFKBIB (By similarity).[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified NKIRAS2 protein (Cat# [TP323273]). The protein was produced from HEK293T cells transfected with NKIRAS2 cDNA clone (Cat# [RC223273]) using MegaTran 2.0 (Cat# [TT210002]).