

## **Product datasheet for TP323277**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Aminoadipate aminotransferase (AADAT) (NM\_016228) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human aminoadipate aminotransferase (AADAT), transcript variant 1,

20 µg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC

>RC223277 representing NM\_016228

or AA Sequence: Red=Cloning site Green=Tags(s)

MNYARFITAASAARNPSPIRTMTDILSRGPKSMISLAGGLPNPNMFPFKTAVITVENGKTIQFGEEMMKR ALQYSPSAGIPELLSWLKQLQIKLHNPPTIHYPPSQGQMDLCVTSGSQQGLCKVFEMIINPGDNVLLDEP AYSGTLQSLHPLGCNIINVASDESGIVPDSLRDILSRWKPEDAKNPQKNTPKFLYTVPNGNNPTGNSLTS ERKKEIYELARKYDFLIIEDDPYYFLQFNKFRVPTFLSMDVDGRVIRADSFSKIISSGLRIGFLTGPKPL IERVILHIQVSTLHPSTFNQLMISQLLHEWGEEGFMAHVDRVIDFYSNQKDAILAAADKWLTGLAEWHVP AAGMFLWIKVKGINDVKELIEEKAVKMGVLMLPGNAFYVDSSAPSPYLRASFSSASPEQMDVAFQVLAQL

**IKESL** 

**TRTRPL**EQKLISEEDLAANDILDYKDDDDK**V** 

Tag: C-Myc/DDK
Predicted MW: 47.2 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Bioactivity:** The specific activity of KATII was determined by measuring the product Kynurenic acid

formation from a conversion of Kynurenine. The reaction was carried out at 37? for 15min in the buffer containing PBS, pH7.4, 2mM a-oxoglutarate, 40µM PLP (pyridoxal 5'-phosphate),

and 0.5mM kynurenine as the substrate.

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.





Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** NP 057312

**Locus ID:** 51166

**UniProt ID:** <u>Q8N5Z0</u>, <u>Q4W5N8</u>

RefSeq Size: 2326 Cytogenetics: 4q33 RefSeq ORF: 1275

Synonyms: KAT2; KATII; KYAT2


**Summary:** This gene encodes a protein that is highly similar to mouse and rat kynurenine

aminotransferase II. The rat protein is a homodimer with two transaminase activities. One activity is the transamination of alpha-aminoadipic acid, a final step in the saccaropine pathway which is the major pathway for L-lysine catabolism. The other activity involves the transamination of kynurenine to produce kynurenine acid, the precursor of kynurenic acid which has neuroprotective properties. Several transcript variants encoding two different


isoforms have been found for this gene. [provided by RefSeq, Nov 2013]

**Protein Pathways:** Lysine biosynthesis, Lysine degradation, Metabolic pathways, Tryptophan metabolism

## **Product images:**







Coomassie blue staining of purified AADAT protein (Cat# TP323277). The protein was produced from HEK293T cells transfected with AADAT cDNA clone (Cat# [RC223277]) using MegaTran 2.0 (Cat# [TT210002]).