

Product datasheet for TP323333


RALY (NM_007367) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RNA binding protein, autoantigenic (hnRNP-associated with lethal yellow homolog (mouse)) (RALY), transcript variant 2
Species:	Human
Expression Host:	HEK293T
Tag:	C-Myc/DDK
Predicted MW:	30.2 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_031393
Locus ID:	22913
RefSeq Size:	1541
Cytogenetics:	20q11.22
RefSeq ORF:	870
Synonyms:	HNRPCL2; P542
Summary:	This gene encodes a member of the heterogeneous nuclear ribonucleoprotein (hnRNP) gene family. This protein may play a role in pre-mRNA splicing and in embryonic development. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]



[View online »](#)

Product images:

Coomassie blue staining of purified RALY protein (Cat# TP323333). The protein was produced from HEK293T cells transfected with RALY cDNA clone (Cat# [RC223333]) using MegaTran 2.0 (Cat# [TT210002]).