


Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001019626</u>
Locus ID:	340526
UniProt ID:	<u>Q5HYW3</u>
RefSeq Size:	4356
Cytogenetics:	Xq13.1
RefSeq ORF:	1707
Synonyms:	6430402L03Rik; MAR5; MART5; RGAG4; SIRH8

Product images:



Coomassie blue staining of purified RTL5 protein (Cat# [TP323344]). The protein was produced from HEK293T cells transfected with RTL5 cDNA clone (Cat# [RC223344]) using MegaTran 2.0 (Cat# [TT210002]).