

Product datasheet for **TP323358L**

NOMO2 (NM_173614) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NODAL modulator 2 (NOMO2), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223358 representing NM_173614 Red =Cloning site Green =Tags(s)

MLVGQGAGLLGPAVVTAAVLLLSGVGPAHGSEDIWVGCGGFVKSDVEINYSLEIKLYTKHGTLYKYQTD
CAPNNGYFMIPLYDKGDFILKIEPLGWSFEPTTVLHVDGVS DICTKGGDINFTGFSVNGKVL SKGQ
PLGPAGVQVSLRNTGTEAKIQSTVTQPGGKFAFFKVLPGDYEILATHPTWALKEASTTVRVTSNANAAS
PLIVAGYNVSGSVRSDGEPMKGVKFLFSSLVTKEDVLGCNVSPVPGFQPQDES LVLYCYTVSREDGSFS
FYSLPSGGYTVIPFYRGERITFDVAPSRLDFTVEHDSLKIEPVFHV MGFSVTGRVLNGPEGDGVPEAVT
LNNQIKVKTKADGSFRLENITTGYTIHAQKEHLYFETVTIKIAPNTPQLADIVATGFSVCGRISIRFP
DTV KQM NKYKVLSSQDKDKSLVTVETDAHGSFCFKAKPGTYKVQVMVPEAETRAGLTLKPQTFPLTVTD
RPVMDVAFVQFLASVSGKVSCLDTCGDLLVTLQSLSRQGEKRLQLSGKVNAMTFTFDNVLPGKYKISIM
HEDWCWKNKSLEVEVLEDDVSAVEFRQTGYMLRCSLSHAITLEFYQDGNGRENVGIYNLSKGVNRFCLSK
PGVYKVT PRSCHRFEQAFYTYDTS SP SILTLTAIRHHVLGTITTDKMMDVTVTIKSSIDSEPALVLGPLK
SVQELRREQQLAEIEARRQEREKNGNEEGEERMTKPPVQEMVDELQGPFSYDFS YWARSGEKITVTPSSK
ELLFYPPSMEAVVSGESCPGKLIEIHGKAGLFLEGQIHPELEGVEIVISEKGASSPLITVFTDDKGAYSV
GPLHSDLEYTVTSQKEGYVLTAVEGTIGDFKAYALAGVSFEIKAEDDQPLPGVLLSLSGGLFRSNLLTQD
NGILTFSNLSPGQYYFKPMMKEFRFEPSSQMIEVQEGQNLKITITGYRTAYSCYGTVSSLNGEPEQGVAM
EAVGQNDCSIYGEDTVTDEEGKFRLRGLLPGCVYHVQLKAEGNDHIERALPHHRVIEVGNNDIDDVNIIV
FRQINQFDLSGNVITSSEYLPTLWVKLYKSENLDNPIQTVSLGQSLFFHFPLL RDGENYVLLDSTLPR
SQYDYILPQVSFTAVGYHKHITLIFNPTRKLPEQDIAQGSYIALPTLLVLLAGYNHDKLIPLLLQLTSR
LQGVGALGQAASDNSGPEDAKRQAKKQKTRRT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	134 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining



[View online »](#)

Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_775885
Locus ID:	283820
UniProt ID:	Q5JPE7
RefSeq Size:	4317
Cytogenetics:	16p12.3
RefSeq ORF:	3666
Synonyms:	Nomo; PM5
Summary:	This gene encodes a protein originally thought to be related to the collagenase gene family. This gene is one of three highly similar genes in a region of duplication located on the p arm of chromosome 16. These three genes encode closely related proteins that may have the same function. The protein encoded by one of these genes has been identified as part of a protein complex that participates in the Nodal signaling pathway during vertebrate development. Mutations in ABCC6, which is located nearby, rather than mutations in this gene are associated with pseudoxanthoma elasticum (PXE). Two transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified NOMO2 protein (Cat# [TP323358]). The protein was produced from HEK293T cells transfected with NOMO2 cDNA clone (Cat# [RC223358]) using MegaTran 2.0 (Cat# [TT210002]).