

Product datasheet for **TP323385**

NOXO1 (NM_172168) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens NADPH oxidase organizer 1 (NOXO1), transcript variant c, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223385 representing NM_172168 Red =Cloning site Green =Tags(s) MAGPRYPVSVQGAALVQIKRLQTFAFSVRWSGDSDTFVRRSWDEFRQLKKTLLKETFVPEAGLLRRSDRVL PKLLGQASLDAPLLGRVGRTSRGLARLQLELYSRLLATAERVARSPITITGFFAPQPLDLEPALPPGSR VILPTPEEQPLSRAAGRLSIHSLEAQLRCLQPFCTQDTRDRPFQAQAQESLDVLLRHPSGWWLVENEDR QTAWFPAPYLEEAAPGQGREGGPSLGSSGPQFCASRAYESSRADELSPAGARVRVLETSDRGWWLCRYG DRAGLLPAVLLRPEGLGALLSGTGFRGGDDPAGEARGFPEPSQATAPPPTVPTRPSPGAIQSRCCTVTRR ALERRRRQGRPRGCVDSVPHPTTEQ SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	41.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_751908</u>



[View online »](#)

Locus ID:	124056
UniProt ID:	Q8NFA2 , A8K832
RefSeq Size:	1147
Cytogenetics:	16p13.3
RefSeq ORF:	1128
Synonyms:	P41NOX; P41NOXA; P41NOXB; P41NOXC; SH3PXD5; SNX28
Summary:	This gene encodes an NADPH oxidase (NOX) organizer, which positively regulates NOX1 and NOX3. The protein contains a PX domain and two SH3 domains. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2012]

Product images:



Coomassie blue staining of purified NOXO1 protein (Cat# TP323385). The protein was produced from HEK293T cells transfected with NOXO1 cDNA clone (Cat# [RC223385]) using MegaTran 2.0 (Cat# [TT210002]).