

## Product datasheet for **TP323398**

### Aldehyde dehydrogenase 10 (ALDH3A2) (NM\_000382) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aldehyde dehydrogenase 3 family, member A2 (ALDH3A2), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223398 representing NM_000382 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MELEVRRVRQAFLSGRSRPLRFRLQQLALRRMVQEREKDLTAIAADLCKSEFNVYSQEVITVLGEIDF  
MLENLPEWVTAKPVKKNVLTMLDEAYIQPQLGVLLIIGAWNYPVLTIQPLIGAIAGNAVIKPELS  
ENTAKILAKLLPQYLDQDLYIVINGGVEETELLKQRFDHIFTGNTAVGKIVMEAAAKHLTPVTLELGG  
KSPCYIDKDCDLDIVCRRITWGYMNCGQTCIAPDYILCEASLQNQIVWKIKETVKEFYGENIKESPDYE  
RIINLRHFKRILSLLEGQKIAFGGETDEATRYIAPTDLTDVDPKTKVMQEEIFGPILPIVPVKNVDEAIN  
FINEREKPLALYVFSHNHKLKRMIDETSSGGVTGNDVIMHFTLNSFPFGGSSGMGAYHGKHSFDTFS  
HQRPCLLKSLKREGANKLRYPPNSQSKVDWGWKFFLLKRFNKEKLGLLLLTLFGIVA AVLKAEYY

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	54.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_000373](#)

Locus ID: 224

UniProt ID: [P51648](#)

RefSeq Size: 3702

Cytogenetics: 17p11.2

RefSeq ORF: 1455


Synonyms: ALDH10; FALDH; SLS

**Summary:** Aldehyde dehydrogenase isozymes are thought to play a major role in the detoxification of aldehydes generated by alcohol metabolism and lipid peroxidation. This gene product catalyzes the oxidation of long-chain aliphatic aldehydes to fatty acid. Mutations in the gene cause Sjogren-Larsson syndrome. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Arginine and proline metabolism, Ascorbate and aldarate metabolism, beta-Alanine metabolism, Butanoate metabolism, Fatty acid metabolism, Glycerolipid metabolism, Glycolysis / Gluconeogenesis, Histidine metabolism, Limonene and pinene degradation, Lysine degradation, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism, Tryptophan metabolism, Valine, leucine and isoleucine degradation

### Product images:



Coomassie blue staining of purified ALDH3A2 protein (Cat# TP323398). The protein was produced from HEK293T cells transfected with ALDH3A2 cDNA clone (Cat# [RC223398]) using MegaTran 2.0 (Cat# [TT210002]).