

## Product datasheet for TP323448L

### NRXN3 (NM\_004796) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human neurexin 3 (NRXN3), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223448 representing NM_004796 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MLGSDDFFVYGGSPSTADLPGSPVSNFMGCLKEVYKNNDIRLELSRLARIADTKMKIYGEVFKCENV ATLDPINFETPEAYISLPKWNTKRMGSISDFRTTEPNGLILFTHGKPQERKDARSQNTKVDFFAVELL DGNLYLLLDMSGTIKVKATQKKANDGEWYHVDIQRDGRSGTISVNSRRTPFASGESEILDLEGDMYLG GLPENRAGLILPTELWTAMLNYGYVCIRDLFIDGRSKNIRQLAEMQNAAGVKSSCSRMSAKQCDSYPCK NNAVCKDGWNRFCICTGTGYWGRTCEREASILSYDGS MYMKIIMPVMHTEAEDVSFRFMSQRAYGLLV ATTSRDSADTLRLELDGGRVKLMVNLDCIRINCNSSKGPETLYAGQKLNDEWHVTVRVRRGKSLKLTVD DDVAEGTMVGDHTRLEFHNIETGIMTEKRYISVVPSSFIGHLQSLMFNGLLYIDLCKNGDIDYCELKARF GLRNIIADPVTFKTKSSYLSLATLQAYTSMHLFFQFKTTPDGFILFNSGDGNDFIAVELVKGYIHVYFD LGNGPNVIKGNDRPLNDNQWHNVITRDNSENTHSLKVDTKVVTVQVINGAKNLDLKGDLVMAGLAQGMYS NLPKLVASRDGFQGLASVDLNGRLPDLINDALHRSGQIERGCEGPSTTCQEDSCANQGVCMQQWEGFTC DCSMTSYSGNQCNPDGATYIFGKSGGLILYTWPANDRPSTRSDRLAVGFSTTVKDGILVRIDSAPGLGDF LQLHIEQGKIGVFNIGTVDISIKEERTPVNDGKYHVVRFRTRNGGNATLQVDNWPVNEHYPTGRQLTIFN TQAQIAIGGKDKGRLFQGLSGLYDGLKVLNMAAENPNIKINGSVRLVGEVPSILGTTQTTSMPPEMS TTVMETTTTMTTTRKNRSTASIQPTSDDLVS SAECSSDDEDFVECEPSTANPTEPGIRRVPGASEVIR ESSSTTGMMVGVAAAAALCILILLYAMYKYRNRDEGSYQVDETRNYISNSAQSNGLMKEKQSSKSGHK KQKNKDREYYV</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	114.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_004787</a>
<b>Locus ID:</b>	9369
<b>UniProt ID:</b>	<a href="#">Q9Y4C0</a> , <a href="#">Q9HDB5</a>
<b>RefSeq Size:</b>	6148
<b>Cytogenetics:</b>	14q24.3-q31.1
<b>RefSeq ORF:</b>	3183
<b>Synonyms:</b>	C14orf60
<b>Summary:</b>	This gene encodes a member of a family of proteins that function in the nervous system as receptors and cell adhesion molecules. Extensive alternative splicing and the use of alternative promoters results in multiple transcript variants and protein isoforms for this gene, but the full-length nature of many of these variants has not been determined. Transcripts that initiate from an upstream promoter encode alpha isoforms, which contain epidermal growth factor-like (EGF-like) sequences and laminin G domains. Transcripts initiating from the downstream promoter encode beta isoforms, which lack EGF-like sequences. Genetic variation at this locus has been associated with a range of behavioral phenotypes, including alcohol dependence and autism spectrum disorder. [provided by RefSeq, Dec 2012]
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Cell adhesion molecules (CAMs)

### Product images:



Coomassie blue staining of purified NRXN3 protein (Cat# [TP323448]). The protein was produced from HEK293T cells transfected with NRXN3 cDNA clone (Cat# [RC223448]) using MegaTran 2.0 (Cat# [TT210002]).