

Product datasheet for TP323654

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Annexin A11 (ANXA11) (NM_145869) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human annexin A11 (ANXA11), transcript variant c, 20 μg

Species: Human
Expression Host: HEK293T

Expression cDNA >RC223654 representing NM_145869
Clone or AA Red=Cloning site Green=Tags(s)

Sequence:

MSYPGYPPPGGYPPAAPGGGPWGGAAYPPPPSMPPIGLDNVATYAGQFNQDYLSGMAANMSGTFGGANM PNLYPGAPGAGYPPVPPGGFGQPPSAQQPVPPYGMYPPPGGNPPSRMPSYPPYPGAPVPGQPMPPPGQQP PGAYPGQPPVTYPGQPPVPLPGQQQPVPSYPGYPGSGTVTPAVPPTQFGSRGTITDAPGFDPLRDAEVLR KAMKGFGTDEQAIIDCLGSRSNKQRQQILLSFKTAYGKDLIKDLKSELSGNFEKTILALMKTPVLFDIYE IKEAIKGVGTDEACLIEILASRSNEHIRELNRAYKAEFKKTLEEAIRSDTSGHFQRLLISLSQGNRDEST NVDMSLAQRDAQELYAAGENRLGTDESKFNAVLCSRSRAHLVAVFNEYQRMTGRDIEKSICREMSGDLEE GMLAVVKCLKNTPAFFAERLNKAMRGAGTKDRTLIRIMVSRSETDLLDIRSEYKRMYGKSLYHDISGDTS

GDYRKILLKICGGND

TRTRPLEQKLISEEDLAANDILDYKDDDDK**V**

Tag: C-Myc/DDK
Predicted MW: 54.2 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some

loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.





RefSeq: NP 665876

Locus ID: 311

UniProt ID: <u>P50995</u>, <u>Q5T0G8</u>


RefSeq Size: 2731 Cytogenetics: 10q22.3 RefSeq ORF: 1515

Synonyms: ALS23; ANX11; CAP-50; CAP50

Summary: This gene encodes a member of the annexin family, a group of calcium-dependent phospholipid-

binding proteins. Annexins have unique N-terminal domains and conserved C-terminal domains, which contain calcium-dependent phospholipid-binding sites. The encoded protein is a 56-kD antigen recognized by sera from patients with various autoimmune diseases. Several transcript variants encoding two different isoforms have been identified. [provided by RefSeq, Dec 2015]

Product images:



Coomassie blue staining of purified ANXA11 protein (Cat# TP323654). The protein was produced from HEK293T cells transfected with ANXA11 cDNA clone (Cat# [RC223654]) using MegaTran 2.0 (Cat# [TT210002]).