

Product datasheet for TP323659M

CCDC83 (NM_173556) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 83 (CCDC83), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223659 representing NM_173556 Red =Cloning site Green =Tags(s)
	<p>MENSGKANKKTDTHDGPPKEIKLPTSEALLDYQCQIKEDAVEQFMFQIKTLRKKKNQKYHERNSRLKEEQIW HIRHLLKELSEEKAEGLPVVTREDVEEAMKEKWKFERDQEKNLDRMRMQISNAEKLFLLEKSEKEYWEEY KNVGSERHAKLITSLQNDINTVKENAEKMSEHYKITLEDTRKKIHKETLLQLDQKKEWATQNAVKLIDKG SYLEIWENDWLKKEVAIHRKEVEELKNAIHELEAENLVLIDQLSNCRLVDLKIPIYPVLHSCPTSNPRHL LLLPLESCLISARRCWRLYLTAAGLEVPEEMSLELPETHIEEKSELQPTEVESRDLMSSEDESTILHL SHENSIEDLQYVKIDKEENSGTEFGDTDMKYLLEYDEKDFKDYVNLGPLGVKLMVESKMPHFQEKI PVKLYKDVRSPESHITYKMMKSFL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	52.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_775827</u>



[View online »](#)

Locus ID: 220047
UniProt ID: [Q8IWF9](#)
RefSeq Size: 2346
Cytogenetics: 11q14.1-q14.2
RefSeq ORF: 1332
Synonyms: CT148; HSD9; KP-CoT-23

Product images:



Coomassie blue staining of purified CCDC83 protein (Cat# [TP323659]). The protein was produced from HEK293T cells transfected with CCDC83 cDNA clone (Cat# [RC223659]) using MegaTran 2.0 (Cat# [TT210002]).