

Product datasheet for TP323701L

PSORS1C2 (NM_014069) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human psoriasis susceptibility 1 candidate 2 (PSORS1C2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223701 representing NM_014069 Red =Cloning site Green =Tags(s)
	MILNWKLLGILVLCLHTRGISGSEGHPSHPPAEDREEAGSPTLPQGPPVPGDPWPGAPPLFEDPPPTRPS RPWRDLPETGVWPPEPPRTDPPQPPRPDDPWAGPQPENPWPPAPEVDNRPQEEDLDPPREEYR
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	14.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_054788
Locus ID:	170680
UniProt ID:	Q9UIG4 , A0A1U9X9A6
RefSeq Size:	1143
Cytogenetics:	6p21.33



[View online »](#)

RefSeq ORF: 408
Synonyms: C6orf17; SPR1
Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified PSORS1C2 protein (Cat# [TP323701]). The protein was produced from HEK293T cells transfected with PSORS1C2 cDNA clone (Cat# [RC223701]) using MegaTran 2.0 (Cat# [TT210002]).