

## Product datasheet for **TP323812M**

### AKAP7 (NM\_138633) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human A kinase (PRKA) anchor protein 7 (AKAP7), transcript variant beta, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223812 representing NM_138633 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MGQLCCFPFSRDEGKISELESSSSAVLQRYSKDIPSWSSGEKNGGEPDDAELVRLSKRLVENAVLKAVQQ YLEETQNKPKPGEGSSVKTEAADQNGNDNENNRK  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	11.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_619539</a>
Locus ID:	9465
UniProt ID:	<a href="#">O43687</a> , <a href="#">Q2TAJ5</a>
RefSeq Size:	2348



[View online »](#)

Cytogenetics: 6q23.2


RefSeq ORF: 312

Synonyms: AKAP15; AKAP18

**Summary:** This gene encodes a member of the A-kinase anchoring protein (AKAP) family, a group of functionally related proteins that bind to a regulatory subunit (RII) of cAMP-dependent protein kinase A (PKA) and target the enzyme to specific subcellular compartments. AKAPs have a common RII-binding domain, but contain different targeting motifs responsible for directing PKA to distinct intracellular locations. Three alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Apr 2011]

**Protein Families:** Druggable Genome

### Product images:



Coomassie blue staining of purified AKAP7 protein (Cat# [TP323812]). The protein was produced from HEK293T cells transfected with AKAP7 cDNA clone (Cat# [RC223812]) using MegaTran 2.0 (Cat# [TT210002]).