

## Product datasheet for **TP323828**

### CPNE1 (NM\_152927) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens copine I (CPNE1), transcript variant 4, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223828 representing NM_152927 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MAHCVTLVQLSISCDHLIDKDIGSKSDPLCVLLQDVGGGSWAELGRTERVRNCSSPEFSKTLQLEYRFET  
VQKLRFGIYDIDNKTPELRDDDGLGGAECSLGQIVSSQVLTLPMLKPGKPAAGRTITVSAQELKDNRW  
TMEVEARNLDKKDFLGKSDPFLEFFRQGDGKWHLVYRSEVIKNNLNPTWKRFVSVQHFCCGPNSTPIQV  
QCSDYDSDGSHDLIGTFHTSLAQLQAVPAEFECIHPEKQKQKKSYSKNSGTIRVKICRVETEYSFLDYVMG  
GCQINFTVGVDFTGSNGDPSSPSLHLYSPTGVNEYLMALWSVGSVWQDYDSDKLFPAFGFGAQVPPDWQ  
VSHEFALNFNPSNPYCAGIQGIVDAYRQALPQVRLYGPTNFAPIINHVARFAAQAHQGTASQYFMLLLL  
TDGAVTDVEATREAVRASNLPM SVIIVGVGGADFEAMEQLDADGGPLHTRSGQAAARDIVQFVPPYRRFQ  
NAPREALAQTVLAEVPTQLVSYFRAQGWAPLKLPPSAKDPAQAPQA

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	58.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP\\_690904](#)

Locus ID: 8904

UniProt ID: [Q99829](#)

RefSeq Size: 2107


Cytogenetics: 20q11.22

RefSeq ORF: 1611

Synonyms: COPN1; CPN1

**Summary:** Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene encodes a calcium-dependent protein that also contains two N-terminal type II C2 domains and an integrin A domain-like sequence in the C-terminus. However, the encoded protein does not contain a predicted signal sequence or transmembrane domains. This protein has a broad tissue distribution and it may function in membrane trafficking. This gene and the gene for RNA binding motif protein 12 overlap at map location 20q11.21. Alternate splicing results in multiple transcript variants encoding different proteins. [provided by RefSeq, Aug 2008]

### Product images:



Coomassie blue staining of purified CPNE1 protein (Cat# TP323828). The protein was produced from HEK293T cells transfected with CPNE1 cDNA clone (Cat# [RC223828]) using MegaTran 2.0 (Cat# [TT210002]).