

Product datasheet for **TP323838L**

BPGM (NM_001724) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human 2,3-bisphosphoglycerate mutase (BPGM), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC223838 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSKYKLIMLRHGEGAWNKENRFCSWVDQKLNSEGMEEARNCGKQLKALNFEFDLVFTSVLNRSIHTAWLI LEELGQEWVPVESSWRLNERHYGALIGLNREQMALNHGEEQVRLWRRSYNVTPPPIEESHPPYYQEIYNDR RYKVC DVPLDQLPRSESLKDVLERLLPYWNERIAPEVLRGKTILISAHGNSRALLKHLEGISDEDIINI TLPTGVPILLELDENLRAVGPHQFLGDQEAIAAIKKVEDQGKVKQAKK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	29.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001715</u>
Locus ID:	669



[View online »](#)


UniProt ID: [P07738](#), [A0A024R782](#)
RefSeq Size: 1800
Cytogenetics: 7q33
RefSeq ORF: 777
Synonyms: DPGM; ECYT8

Summary: 2,3-diphosphoglycerate (2,3-DPG) is a small molecule found at high concentrations in red blood cells where it binds to and decreases the oxygen affinity of hemoglobin. This gene encodes a multifunctional enzyme that catalyzes 2,3-DPG synthesis via its synthetase activity, and 2,3-DPG degradation via its phosphatase activity. The enzyme also has phosphoglycerate phosphomutase activity. Deficiency of this enzyme increases the affinity of cells for oxygen. Mutations in this gene result in hemolytic anemia. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Sep 2009]

Protein Families: Druggable Genome

Protein Pathways: Glycolysis / Gluconeogenesis, Metabolic pathways

Product images:



Coomassie blue staining of purified BPGM protein (Cat# [TP323838]). The protein was produced from HEK293T cells transfected with BPGM cDNA clone (Cat# [RC223838]) using MegaTran 2.0 (Cat# [TT210002]).