

## Product datasheet for **TP323980M**

### **XRN2 (NM\_012255) Human Recombinant Protein**

#### **Product data:**

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human 5'-3' exoribonuclease 2 (XRN2), 100 µg

**Species:** Human

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >RC223980 representing NM\_012255

**Red**=Cloning site **Green**=Tags(s)

MGVPAFFRWLSRKYPPIVNCVEEKPKCEKNGVKIPVDASKPNPNDVEFDNLYLDMNGIIHPCTHPEDKPA  
PKNEDEMMVAIFEYIDRLFSIVRPRRLLYMAIDGVAPRAKMNQQRSSRRFRASKEGMEAAVEKQVRVEEIL  
AKGGFLPPEEIKERFDSNCITPGTEFMDNLAKCLRYIADRLNNDPGWKNLTVILSDASAPGEGEHKIMD  
YIRRQRAQPNHDPNTHHCLCGADADLIMLGLATHEPNFTIREEFKPNKPKCGLCNQFGHEVKDCEGLP  
REKKGKHDELADSLPCAEGEFIFLRLNVLREYLERELTMSALPFTFDVERSIDDWVFMCFVGNDFLPHL  
PSLEIRENAIDRLVNIYKNVHKTGGYLTESGYVNLQRVQMIMLAVGEVEDSIFKKRKDDSDSFRRRQKE  
KRKRMKRDQPAFTPSGILTPHALGSRNSPGSQVASNPRQAAYEMRMQNNSSPSISPNTSFTSDGSPSPLG  
GIKKAEDSDSEPEPEDNVRLWEAGWKQRYKKNKFDVDAADEKFRKVVQSYVEGLCWVLRYYYQGCASW  
KWYYPFHYPFASDFEGIADMPSDFEKGTKPFKPLEQLMGVFPAAASGNFLPPSWRKLMSDPDSSIIDFYP  
EDFAIDLNGKKYAWQGVALLPFVDERRLRAALEEVYDLPTEETRRNSLGGDVLVFGKHHPLHDFILELY  
QTGSTPEVPELCHGIQGFSLDEEAILPDQIVCSPVPMRLDLTQNTVVSINFKDPQFAEDYIFKAVM  
LPGARKPAAVLKPSPDWEKSSNGRQWKPQLGFNRDRRPVHLDQAAFRTLGHVMPRGSGTGIYSNAAPPVPT  
YQGNLYRPLLRGQAQIPKLMNMRPQDSWRGPPPLFQQQRFDRGVGAEP LLPWNRMLQTQNAAFQPNQYQ  
MLAGPGGYPPRRDDRGRQGYPREGRKYPLPPPSGRYNWN

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK

**Predicted MW:** 108.4 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol


**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_036387</a>
<b>Locus ID:</b>	22803
<b>UniProt ID:</b>	<a href="#">Q9H0D6</a>
<b>RefSeq Size:</b>	3449
<b>Cytogenetics:</b>	20p11.22
<b>RefSeq ORF:</b>	2850
<b>Summary:</b>	This gene encodes a 5'-3' exonuclease that promotes transcription termination at cotranscriptional cleavage sites. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015]
<b>Protein Pathways:</b>	RNA degradation

### Product images:



Coomassie blue staining of purified XRN2 protein (Cat# [TP323980]). The protein was produced from HEK293T cells transfected with XRN2 cDNA clone (Cat# [RC223980]) using MegaTran 2.0 (Cat# [TT210002]).