

Product datasheet for **TP324011M**

PEDS1 (NM_199129) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transmembrane protein 189 (TMEM189), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224011 representing NM_199129 Red =Cloning site Green =Tags(s) MAGAEDWPGQQLELDEDEASCCRWGAQHAGARELAALYSPGKRLQEWCSVILCFSLIAHNLVHLLLLAR W EDTPLVILGVVAGALIADFLSGLVHWGADTWGSVELPIVGKAFIRPFREHHIDPTAIRHDFIETNGDNC LVTLLPLLNMAYKFRTHSPEALEQLYPWECFVFCLIFGTFTNQIHKWSHTYFGLPRWVTLLQDWHVILP RKHHRIHHVSPHETYFCITTGWLNYPLEKIGFWRRLEDLIQGLTGEKPRADDMKWAQKIK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	31 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_954580</u>
Locus ID:	387521


[View online »](#)

UniProt ID:	<u>A5PLL7</u>
RefSeq Size:	2187
Cytogenetics:	20q13.13
RefSeq ORF:	810
Synonyms:	CarF; KUA; TMEM189
Summary:	Co-transcription of this gene and the neighboring downstream gene (ubiquitin-conjugating enzyme E2 variant 1) generates a rare read-through transcript, which encodes a fusion protein comprised of sequence sharing identity with each individual gene product. The protein encoded by this individual gene lacks a UEV1 domain but includes three transmembrane regions. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2009]
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified TMEM189 protein (Cat# [TP324011]). The protein was produced from HEK293T cells transfected with TMEM189 cDNA clone (Cat# [RC224011]) using MegaTran 2.0 (Cat# [TT210002]).