

Product datasheet for **TP324090L**

QK1 (QKI) (NM_206853) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human quaking homolog, KH domain RNA binding (mouse) (QKI), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224090 representing NM_206853 Red =Cloning site Green =Tags(s)

MVGEMETKEKPKPTPDYLMQLMNDKKLMSSLPNFCGIFNHLERLLDDEISRVRKDMYNDTLNGSTEKRSA
ELPDAVGPIVQLQEKLIVPVKEYPDFNFVGRILGPRGLTAKQLEAETGCKIMVRGKGSMDKKKEEQNRG
KPNWEHLNEDLHVLITVEDAQNRAEIKLRAVEEVKLLVPAEGEDSLKMKQLMELAILNGTYRDANIK
SPALAFSLAATAQAAPRIITGPAPVLPAAALRTPPTAGPTIMPLIRQIQTAVMPNGTPHPTAAIVPPGPE
AGLIYTPYEYPYTLAPATSILEYPIEPSGVLGMAFPTKG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	35 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_996735
Locus ID:	9444



[View online »](#)

UniProt ID:	Q96PU8 , Q8WY44
RefSeq Size:	5693
Cytogenetics:	6q26
RefSeq ORF:	957
Synonyms:	Hqk; hqkl; QK; QK1; QK3
Summary:	The protein encoded by this gene is an RNA-binding protein that regulates pre-mRNA splicing, export of mRNAs from the nucleus, protein translation, and mRNA stability. The encoded protein is involved in myelination and oligodendrocyte differentiation and may play a role in schizophrenia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]

Product images:



Coomassie blue staining of purified QKI protein (Cat# [TP324090]). The protein was produced from HEK293T cells transfected with QKI cDNA clone (Cat# [RC224090]) using MegaTran 2.0 (Cat# [TT210002]).