

Product datasheet for **TP324133L**

PDE9A (NM_001001568) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens phosphodiesterase 9A (PDE9A), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC224133 representing NM_001001568
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MDAFRSTPYKVRPVAIKQLSEREELIQSVLAQVAEQFSRAFKINELKAEVANHLAVLEKRVELEGLKVVE
IEKCKSDIKKMREELAARSSRTNCPCKYSFLDNHKKLTPRRDVPTYPKYLLSPETIEALRKPTFDVWLWE
PNEMLSCLEHMYHDLGLVRDFSINPVTLRRWLFVHDNYRNNPFHNFRHCFVAQMMYSMVWLCSLQKEF
SQTDLILMTAAICHDLDPGYNNNTYQINARTELVRYNDISPLENHHCAVAFQILAEPECNIFSNIPDP
GFKQIRQGMITLILATDMARHAEIMDSFKEKMENFDYSNEEHMTLLKMILIKCCDISNEVRPMEVAEPWW
DCLLEEFMQSDREKSEGLPVAPFMDRDKVTKATAQIGFIKVLIPMFETVTKLFPMVEEIMLQPLWESR
DRYEELKRIDDAMKELQKKTDSLTSGATEKSRERSRDVKNSEGDC

TRRLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	54.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP_001001568](#)

Locus ID: 5152

UniProt ID: [O76083](#)

RefSeq Size: 1774

Cytogenetics: 21q22.3

RefSeq ORF: 1398


Synonyms: HSPDE9A2

Summary: The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Progesterone-mediated oocyte maturation, Purine metabolism

Product images:



Coomassie blue staining of purified PDE9A protein (Cat# [TP324133]). The protein was produced from HEK293T cells transfected with PDE9A cDNA clone (Cat# [RC224133]) using MegaTran 2.0 (Cat# [TT210002]).