

## Product datasheet for TP324199

### C19orf54 (NM\_198476) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 19 open reading frame 54 (C19orf54), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224199 representing NM_198476 Red=Cloning site Green=Tags(s)
	<p>MTSPCSPPLKPPISPPKTPVPQASSIPSPPLPSPPLDFSALPSPPPWSQQTVPVPPPLPLPPPPAATGPAPR            HVFGLEKSQLLKEAFDKAGPVPKGREVDKRLKLLKHKDRFRGDLRWILFCADLPSLIQEGPQCGLVALWMA            GTLLSPPSGVPLERLIRVATERGYTAQGEMFSVADMGRLAQEVLCQAKLLSGGLGGPNRDLVLQHLVTG            HPLLIPYDEDFNHEPCQRKKGHKAHWAVSAGVLLGVRAPSLGYTEDPELPLGHPVLGTPCQPPSLPEEG            SPGAVYLLSKQGKSWHYQLWDYDQVRESNLQLTDFSPSRATDGRVYVVPVGGVRAGLCGQALLLTPQDCS            H</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	37.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_940878</a>
Locus ID:	284325



[View online »](#)

UniProt ID: [Q5BKX5](#)  
RefSeq Size: 2769  
Cytogenetics: 19q13.2  
RefSeq ORF: 1053

### Product images:



Coomassie blue staining of purified C19orf54 protein (Cat# TP324199). The protein was produced from HEK293T cells transfected with C19orf54 cDNA clone (Cat# [RC224199]) using MegaTran 2.0 (Cat# [TT210002]).