

## Product datasheet for TP324235

### PAOX (NM\_207128) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human polyamine oxidase (exo-N4-amino) (PAOX), transcript variant 5, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224235 representing NM_207128 Red=Cloning site Green=Tags(s)

MESTGSVGEAPGGPRVLVGGGIAGLGAAQRLCGHSAPHLRVLEATARAGGRIRSERCFGGWVEVGAHW  
IHGPSRGNPVFQLAAEYGLLGEKELSQENQLVETGGHVGLPSVSYASSGASVSLQVAEMATLFYGLIDQ  
TREFLHAAETPVPSVGEYLLKKEIGQHVAGWTEDETRKLLAVLNSFFNLECCVSGTHSMDLVALAPFGE  
YTVLPGLDCTFSKGYQGLTNCMMAALPEDTVVFEKPKVTIHWNGSFQEAAFPGETFPVSVECEDGDRFPA  
HHVIVTVPLGFLREHLDTFFDPPLPAEKAEAIRKIGFGTNNKIFLEFEFPWEPDCQLIQLVWEDTSPL  
DAAPQLQDAWFRKLIGFVVLPAFASVHVLGFIAGLESEFMETLSDEEVLLCLTQVLRRTAPDPVCGGS  
HTSHVLLHDARGSAVGMEGGRPPPQSVGPAGAAAQAQALAGPSLLCSTRVGGRLGPSFLLTDFSLA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	51.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_997011](#)

Locus ID: 196743

UniProt ID: [Q6QHF9](#)

RefSeq Size: 1720

Cytogenetics: 10q26.3


RefSeq ORF: 1458

Synonyms: PAO

**Summary:** Flavoenzyme which catalyzes the oxidation of N(1)-acetylspermine to spermidine and is thus involved in the polyamine back-conversion (PubMed:12477380). Can also oxidize N(1)-acetylspermidine to putrescine. Substrate specificity: N(1)-acetylspermine = N(1)-acetylspermidine > N(1),N(12)-diacetylspermine >> spermine. Does not oxidize spermidine. Plays an important role in the regulation of polyamine intracellular concentration and has the potential to act as a determinant of cellular sensitivity to the antitumor polyamine analogs (PubMed:12477380).[UniProtKB/Swiss-Prot Function]

**Protein Families:** Druggable Genome

### Product images:



Coomassie blue staining of purified PAOX protein (Cat# TP324235). The protein was produced from HEK293T cells transfected with PAOX cDNA clone (Cat# [RC224235]) using MegaTran 2.0 (Cat# [TT210002]).