

Product datasheet for **TP324309L**

C14orf181 (NM_207442) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 14 open reading frame 181 (C14orf181), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224309 representing NM_207442 Red =Cloning site Green =Tags(s) MNTYTRSASFPTKFTFPLAQNRTSHPTTTSKLYLGRPSRRDALQANGEGPEVGPTEEPAHRLPVGS GSAGEAAASLGTARVRSGRAECTAKLAGDGHARLASAFLHRRALPPPLGTCPRLGEGGSAPALGPSPR RGRRVSEGLRDSLPLPLGSLTLGKSPK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_997325
Locus ID:	400223
RefSeq Size:	1727
Cytogenetics:	14q24.1




[View online »](#)

RefSeq ORF: 498

Synonyms: FLJ39779

Summary: DISCONTINUED: This record has been withdrawn by NCBI and HGNC. This record was defined by AK097098.1, which is not sufficient evidence to define a distinct locus.

Product images:



Coomassie blue staining of purified C14orf181 protein (Cat# [TP324309]). The protein was produced from HEK293T cells transfected with C14orf181 cDNA clone (Cat# [RC224309]) using MegaTran 2.0 (Cat# [TT210002]).