

## Product datasheet for TP324314

### SSX4 (NM\_005636) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human synovial sarcoma, X breakpoint 4 (SSX4), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC224314 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MNGDDAFARRPRDDAQISEKLRKAFDDIAKYFSKKEWEKMKSSSEKIVVYVMKLNIEVMTKLGFKVTLPPF  
MRSKRAADFHGNDFGNDRNHRNQVERPQMTFGSLQRIFPKIMPKKPAEEENGLKEVPEASGPQNDGKQLC  
PPGNPSTLEKINKTSGPKRGKHAWTHRLRERKQLVVYEEISDPEDDE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	21.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_005627</a>
Locus ID:	6759
UniProt ID:	<a href="#">O60224</a>



[View online >](#)

RefSeq Size: 1250

Cytogenetics: Xp11.23


RefSeq ORF: 564

Synonyms: CT5.4

**Summary:** The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneously humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy. SSX1, SSX2 and SSX4 genes have been involved in the t(X;18) translocation characteristically found in all synovial sarcomas. This translocation results in the fusion of the synovial sarcoma translocation gene on chromosome 18 to one of the SSX genes on chromosome X. Chromosome Xp11 contains a segmental duplication resulting in two identical copies of synovial sarcoma, X breakpoint 4, SSX4 and SSX4B, in tail-to-tail orientation. This gene, SSX4, represents the more telomeric copy. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:** Transcription Factors

### Product images:



Coomassie blue staining of purified SSX4 protein (Cat# TP324314). The protein was produced from HEK293T cells transfected with SSX4 cDNA clone (Cat# [RC224314]) using MegaTran 2.0 (Cat# [TT210002]).