

Product datasheet for **TP324407**

N42L1 (N4BP2L1) (NM_052818) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NEDD4 binding protein 2-like 1 (N4BP2L1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224407 representing NM_052818 Red =Cloning site Green =Tags(s)

MEDSFLQSFGRSLQPQQQQRQRPPRPPRGTPPPRRHSFRKHLYLLRGLPGSGKTTLARQLQHDFPRAL
IFSTDDFFFREDGAYEFNPDFLEEAHEWNQKRARKAMRNGISPIIIDNTNLHAWEMKPYAVMALENNYEV
IFREPDTRWKFNVQELARRNIHGVSREKIHRMKERYEHDVTFHSVLHAEKPSRMNRNQDRNNALPSNNAR
YWNSYTEFPNRRAHGGFTNESSYHRRGGCHHG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	28.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_438169
Locus ID:	90634



[View online »](#)

UniProt ID: [Q5TBK1](#)
RefSeq Size: 3051
Cytogenetics: 13q13.1
RefSeq ORF: 729
Synonyms: CG018

Product images:

Coomassie blue staining of purified N4BP2L1 protein (Cat# TP324407). The protein was produced from HEK293T cells transfected with N4BP2L1 cDNA clone (Cat# [RC224407]) using MegaTran 2.0 (Cat# [TT210002]).