

## Product datasheet for TP324474

### ARF1 (NM\_001024226) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ADP-ribosylation factor 1 (ARF1), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224474 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MGNIFANLFKGLFGKKEMRILMVGLDAAGKTTILYKLLGGEIVTTIPTIGFNVETVEYKNISFTWWDVGG QDKIRPLWRHYFQNTQGLIFVWDSNDRERVNEAREELMRMLAEDELRDVLLVFANKQDLPNAMNAAEI T DKLGLHSLRHRNWYIQATCATSGDGLYEGLDWLSNQLRNQK  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	20.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001019397</a></u>
Locus ID:	375
UniProt ID:	<u><a href="#">P84077</a></u>



[View online »](#)

RefSeq Size: 1986

Cytogenetics: 1q42.13


RefSeq ORF: 543

Synonyms: PVNH8

**Summary:** ADP-ribosylation factor 1 (ARF1) is a member of the human ARF gene family. The family members encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking as activators of phospholipase D. The gene products, including 6 ARF proteins and 11 ARF-like proteins, constitute a family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2 and ARF3), class II (ARF4 and ARF5) and class III (ARF6), and members of each class share a common gene organization. The ARF1 protein is localized to the Golgi apparatus and has a central role in intra-Golgi transport. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]

**Protein Pathways:** Vibrio cholerae infection

### Product images:



Coomassie blue staining of purified ARF1 protein (Cat# TP324474). The protein was produced from HEK293T cells transfected with ARF1 cDNA clone (Cat# [RC224474]) using MegaTran 2.0 (Cat# [TT210002]).