

Product datasheet for TP324481L

OGT (NM_181673) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human O-linked N-acetylglucosamine (GlcNAc) transferase (UDP-N-acetylglucosamine:polypeptide-N-acetylglucosaminyl transferase) (OGT), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224481 representing NM_181673 Red=Cloning site Green=Tags(s)

MASSVGNVADSTGLAELAHREYQAGDFEAAERHCMQLWRQEPDNTGVLLLLSSIHFCRRLDRSAHFSTL
AIKQNPLLAAYSNLGNVYKERGQLQEAEHYRHALRLKPDFIDGYINLAAALVAAGDMEGAVQAYVSAL
QYNPDLYCVRSDLGNLLKALGRLEEAKACYLKAJETQPNFAVAWSNLGCVFNAQGEIWLAIHHEKAVTL
DPNFLDAYINLGNVLKEARIFDRAVAAYLRALSLSPNHAWVHGNLACVYEQGLIDLDAIDTYRRAIELQP
HFPDAYCNLANALKEKGSVAEAEDCYNTALRLCPTHADSLNLANIKREQGNIEEAVRLYRKALEVFPEF
AAHSNLA SVLQQGKLQEALMHYKEAIRISPTFADAYS NMGNLTKEMQDVQGALQCYTRAIQINPAFAD
AHSNLA SIHKDSGNIPEAIASYRTALKLKPDPDAYCNLAHCLQIVCDWTDYDERMKKLVSIVADQLEKN
RLPSVHPHHSMLYPLSHGFRKAIAERHGNLCLDKINVLHKPPYEHPKDLKLSDGRLRVGYVSSDFGNHPT
SHLMQSIPGMHNPDKFEVFCYALSPDDGTNFRVKVMAEANHFIDLSQIPCNGKAADRIHQDGIHILVNMN
GYTKGARNELFALRPAIQAMWLGYPGTSGALFMDYIITDQETSPAEEAEQYSEKLAYMPHTFFIGDHAN
MFPHLKKKAVIDFKSNGHIYDNRIVLNGIDLKAFDLSLPVKIVKMKCPDGGDNADSSNTALNMPVPMN
TIAEAVIEMINRGQIQITINGFSISNGLATTQINNKAATGEEVPTIIVTTRSQYGLPEDAIVYCNFNQL
YKIDPSTLQMWANILKRVPN SVLWLLRFAVGEPNIQQAQNMGLPQNRIIFSPVAPKEEHVRRGQLADV
CLDTPLCNGHTTGMDVLWAGTPMVTMPGETLASRVAASQLTCLGCLELIAKNRQEYEDIAVKLGTDLEYL
KKVRGKVVWKQRISPLFNTKQYTMELERLYLQMWHEYAAGNKPDMIKPVEVTESA

SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	115.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_858059
Locus ID:	8473
UniProt ID:	O15294
RefSeq Size:	5445
Cytogenetics:	Xq13.1
RefSeq ORF:	3108
Synonyms:	HINCUT-1; HRNT1; MRX106; O-GLCNAC; OGT1
Summary:	This gene encodes a glycosyltransferase that catalyzes the addition of a single N-acetylglucosamine in O-glycosidic linkage to serine or threonine residues. Since both phosphorylation and glycosylation compete for similar serine or threonine residues, the two processes may compete for sites, or they may alter the substrate specificity of nearby sites by steric or electrostatic effects. The protein contains multiple tetratricopeptide repeats that are required for optimal recognition of substrates. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Oct 2009]
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, O-Glycan biosynthesis

Product images:



Coomassie blue staining of purified OGT protein (Cat# [TP324481]). The protein was produced from HEK293T cells transfected with OGT cDNA clone (Cat# [RC224481]) using MegaTran 2.0 (Cat# [TT210002]).