

Product datasheet for **TP324570**

SCGB1C2 (NM_001097610) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human secretoglobin, family 1C, member 1-like (LOC653486), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224570 representing NM_001097610 Red =Cloning site Green =Tags(s) MKGSRALLLVALTLFCICRMATGEDNDEFFMDFLQTLVGTPEELYEGTLGKYNVNEDAKAAMTELKSCR DGLQPMHKAELVKLLVQVLGSQDGA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	10.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001091079
Locus ID:	653486
UniProt ID:	Q8TD33 , P0DMR2
RefSeq Size:	436
Cytogenetics:	17p13.3




[View online »](#)

RefSeq ORF: 285

Synonyms: SCGB1C1

Product images:



Coomassie blue staining of purified SCGB1C2 protein (Cat# TP324570). The protein was produced from HEK293T cells transfected with SCGB1C2 cDNA clone (Cat# [RC224570]) using MegaTran 2.0 (Cat# [TT210002]).