

Product datasheet for TP324816

WNT8A (NM_058244) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human wingless-type MMTV integration site family, member 8A (WNT8A), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224816 representing NM_058244 Red =Cloning site Green =Tags(s)

MGNLFMLWAALGICCAAFSASAWSVNNFLITGPKAYLTYTTSVALGAQSGIEECKFQFAWERWNCPENAL
QLSTHNRLRSATRETSFIHAISSAGVMYIITKNCSMGDFENCDCDGSNNGKTGGHGWIWGGSDNVEFG
E
RISKLFVDSLEKGDARALMNLHNNRAGRLAVRATMKRTCKCHGISGSCSIQTCWLQLAEFREMGDYLKA
KYDQALKIEMDKRQLRAGNSAEGHWVPAEAFPSAEALIFLEESPDYCTCNSSLGIYGTEGRECLQNSH
NTSRWERRSCGRLCTEGLQVEERKTEVISSCNCKFQWCCTVKCDQCRHVSVKYYCARSPGSAQSLGKGS
A

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	38.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP_490645](#)

Locus ID: 7478

UniProt ID: [Q9H1J5](#)

RefSeq Size: 1899

Cytogenetics: 5q31.2

RefSeq ORF: 1053


Synonyms: WNT8D

Summary: The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family, and may be implicated in development of early embryos as well as germ cell tumors. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2014]

Protein Families: Cancer stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Secreted Protein, Stem cell relevant signaling - Wnt Signaling pathway

Protein Pathways: Basal cell carcinoma, Hedgehog signaling pathway, Melanogenesis, Pathways in cancer, Wnt signaling pathway

Product images:



Coomassie blue staining of purified WNT8A protein (Cat# TP324816). The protein was produced from HEK293T cells transfected with WNT8A cDNA clone (Cat# [RC224816]) using MegaTran 2.0 (Cat# [TT210002]).