

Product datasheet for TP324846

ZNF238 (ZBTB18) (NM_205768) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 238 (ZNF238), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224846 protein sequence Red=Cloning site Green=Tags(s)

MCPKGYEDSMFEPPDHSRHLLQCLSEQRHQGFCDCTVLVGDAQFRAHRAVLASCSMYFHIFYKDQLDKRD
IVHLNSDIVTAPAFALLLEFMYEGKLQFKDLPIEDVLAASYLHMYDIVKVCKKLEKATTEADSTKKE
EDASSCDKVESLSDGSSHIAGDLPSEDEGEDEKLNILPSKRDAAEPGNMWMRLPDSAGIPQAGGEE
EPHATAAGKTVASPCSSTESLSQRSVTSVRDSADVDCVLDLSVKSSLSGVENLNSSYFSSQDVLRSNLVQ
VKVEKEASCDSDVGTNDYDMEHSTVKESVSTNNRVQYEPHLLAPLREDSVLRDREDKASDDEMTPPE
SERVQVEGGMESLLPYVSNILSPAGQIFMCPKCNKVFSPHILQIHLSTHFREQDGIRSKPAADVNVPT
CSLCGKTFSCMYTLKRHERTHSGEKPYTCTQCGKSFQYSHNLSRHAVVHTREKPHACKWCERRFTQSGDL
YRHIRKHFCELVNSLSVKSEALSPTVRDWTLEDSSQELWK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	59.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_991331](#)

Locus ID: 10472

UniProt ID: [Q99592](#)

RefSeq Size: 3922

Cytogenetics: 1q44


RefSeq ORF: 1593

Synonyms: C2H2-171; MRD22; RP58; TAZ-1; ZNF238

Summary: This gene encodes a C2H2-type zinc finger protein which acts a transcriptional repressor of genes involved in neuronal development. The encoded protein recognizes a specific sequence motif and recruits components of chromatin to target genes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified ZBTB18 protein (Cat# TP324846). The protein was produced from HEK293T cells transfected with ZBTB18 cDNA clone (Cat# [RC224846]) using MegaTran 2.0 (Cat# [TT210002]).