

Product datasheet for TP324944M

RANBP3 (NM_007320) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RAN binding protein 3 (RANBP3), transcript variant RANBP3-b, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC224944 representing NM_007320 Red =Cloning site Green =Tags(s)

MADLANEEKPAIAPPVVFVQKDKGQKRSAGGSSPEGGEDSDREDGNYCPPVKRERTSSLTQFPSPQSEER
SSGFRLKPPTLIHQAPSAGLPSQKPKEQQRSVLRPAVLQAPQPKALSQTVPSSGTNGVSLPADCTGAVP
AASPDAAWRSPSEAADEVCALEEKPEQKNESSNASEEEACEKDPATQQAQAFVFGQNLDRVVKLINESVD
EADMENAGHPSADTPTATNYFLQYISSLENSTNSADASSNKVFGQNMSEVLSPPKLNEVSSDANREN
AAAESGSESPFQEATPEKESLAESAAAYTKATARKCLLEKVEVITGEEAESNLVQMCKLFVFDKTSQSW
VERGRLLRLNDMASTDDGTLQSRLVMRTQGSLRLILNTKLWAQMIDKASEKSIRITAMDTEDQGKVF
LISASKDTGQLYAALHHRILALRSRVEQEQAEMPAPEPGAAPSNEEDSDDDDDVLAPSGATAAGAGDE
GDGQTTGST

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	53.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_015559
Locus ID:	8498
UniProt ID:	Q9H6Z4 , Q9H6Z4-3
RefSeq Size:	1868
Cytogenetics:	19p13.3
RefSeq ORF:	1497
Summary:	This gene encodes a protein with a RanBD1 domain that is found in both the nucleus and cytoplasm. This protein plays a role in nuclear export as part of a heteromeric complex. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified RANBP3 protein (Cat# [TP324944]). The protein was produced from HEK293T cells transfected with RANBP3 cDNA clone (Cat# [RC224944]) using MegaTran 2.0 (Cat# [TT210002]).