

## Product datasheet for **TP324956L**

### Insulin Receptor R (INSRR) (NM\_014215) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human insulin receptor-related receptor (INSRR), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC224956 representing NM_014215
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MAVPSLWPWGACLPVIFLSLGFGLDVEVCPSLDIRSEVAELRQLENCVVEGHLQILLMFTATGEDFRG  
LSFPRLTQVTDYLLLFRVYGLESRLDLPNLAIVIRGTRFLGYALVIFEMPHLRDVALPALGAVLRGAVR  
VEKNQELCHLSTIDWGLLQPAPGANHIVGNKLGEECADVCPGVLGAAGEPCAKTTFSGHTDYRCWTSSHC  
QRVCPCPHGMACTARGECCHECLGGCSQPEDPRACVACRHLVFQGAQLWACPPGTQYQESWRCVTAERC  
ASLSVPGRASTFGIHQGSCLAQCPSGFTRNSSSIFCHKCEGLCPKECKVGTKTIDSIQAAQDLVGCTHV  
EGSLILNLRQGYNLEPQLQHSLGLVETITGFLKIKHSFALVSLGFFKNLKLIRGDAMVDGNYTLYVLDNQ  
NLQQLGSWVAAGLTIPVGKIYFAFNPRLCLEHIYRLEEVTGTRGRQNKAEINPRTNGDRAACQTRTLRFV  
SNVTEADRILLRWERYEPEARLDSFIVYKESPFQNAHEVGPDACGTQSWNLLDVELPLSRTQEPGV  
TLASLKPWTQYAVFVRAITLTTEEDSPHQGAQSPIVYLRLTLPAAPTVPQDVISTSNSSSHLLVRWKPTQ  
RNGNLTYLVLWQRLAEDGDLYLNDYCHRLRLPTSNNDRFDGEDGDPEAEMESDCCPCQHPPPGQVLP  
PLEAQEASFQKKFENFLHNAITIPISPKVTSINKSPQRDSGRHRRRAAGPLRLGGNSSDFEIQEDKVPRE  
RAVLSGLRHFTEYRIDIHACNHAHTVGCSAATFVFARTMPHREADGIPGKVAWEASSKNSVLLRWLEPP  
DPNGLILKYEIKYRRLGEEATVLCVSRLRYAKFGGVHLALLPPGNYSARVRATSLAGNGSWTDSVAFYIL  
GPEEDAGGLHVLLTATPVGLTLLIVLAALGFFYGKKRNRTLYASVNPEYFSASDMYVPDEWEVPREQIS  
IIRELGQGSFGMVYEGLARGLEAGEESTPVALKTVNELASPRECIEFLKEASVMKAFKCHHVRLGVS  
QGQPTLVIMELMTRGDLKSHLRLRPEAENNPGLPQPALGEMIQMAGEIADGMAYLAANKFVHRDLAARN  
CMVSQDFTVKIGDFGMTRDVYETDYRKGKGLLPVRWMAPELKDGIPTTHSDVWSFGVWLWEIVTLAE  
QPYQGLSNEQVLKFVMDGGVLEELEGCPQLQLQELMSRCWQPNPRLRPSFTHILDSIQEELRPSFRLLSFY  
YSPECRGARGSLPTTDAEPDSSPTPRDCSPQNGGPGH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	140.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method



[View online »](#)

<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_055030</a>
<b>Locus ID:</b>	3645
<b>UniProt ID:</b>	<a href="#">P14616</a>
<b>RefSeq Size:</b>	4148
<b>Cytogenetics:</b>	1q23.1
<b>RefSeq ORF:</b>	3891
<b>Synonyms:</b>	IRR
<b>Summary:</b>	Receptor with tyrosine-protein kinase activity. Functions as a pH sensing receptor which is activated by increased extracellular pH. Activates an intracellular signaling pathway that involves IRS1 and AKT1/PKB.[UniProtKB/Swiss-Prot Function]
<b>Protein Families:</b>	Druggable Genome, Protein Kinase, Transmembrane
<b>Protein Pathways:</b>	Prostate cancer, Regulation of actin cytoskeleton

### Product images:



Coomassie blue staining of purified INSRR protein (Cat# [TP324956]). The protein was produced from HEK293T cells transfected with INSRR cDNA clone (Cat# [RC224956]) using MegaTran 2.0 (Cat# [TT210002]).