

Product datasheet for TP325059L

YPEL5 (NM_001127399) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens yippee-like 5 (Drosophila) (YPEL5), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225059 protein sequence Red =Cloning site Green =Tags(s)
	MGRIFLDHIGGTRLFSCANCDTILNRSSELISTRFTGATGRAFLFNKWNLQYSEVQDRVMLTGRHMVRD VSCKNCNSKLGWIYEFATEDSQRYKEGRVILERALVRESEGFEEHVPSDNS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001120871</u>
Locus ID:	51646
UniProt ID:	<u>P62699</u>
RefSeq Size:	2342



[View online »](#)


Cytogenetics: 2p23.1

RefSeq ORF: 363

Synonyms: CGI-127

Summary: Component of the CTLH E3 ubiquitin-protein ligase complex that selectively accepts ubiquitin from UBE2H and mediates ubiquitination and subsequent proteasomal degradation of the transcription factor HBP1 (PubMed:29911972). Required for normal cell proliferation (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified YPEL5 protein (Cat# [TP325059]). The protein was produced from HEK293T cells transfected with YPEL5 cDNA clone (Cat# [RC225059]) using MegaTran 2.0 (Cat# [TT210002]).