

Product datasheet for **TP325113**

FAM136B (NM_001012983) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 136, member B (FAM136B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225113 representing NM_001012983 Red=Cloning site Green=Tags(s)

MVELQQLRVQEVVDSMVKSLERENIWKTQGLMFWCSASCCEDSQAFTQQVHQIECCPVPLAQVQALVTS
ELEKFQDHLARCTMHCNDKAKDSIDAGSKELQVKQQLDGCVTCKVDDHHLISTMTKKMKEAVLSIGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	15.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001013001
Locus ID:	387071
Cytogenetics:	6p25.2
RefSeq ORF:	414



[View online »](#)

Synonyms: C6orf87; dj40E16.3

Product images:



Coomassie blue staining of purified FAM136B protein (Cat# TP325113). The protein was produced from HEK293T cells transfected with FAM136B cDNA clone (Cat# [RC225113]) using MegaTran 2.0 (Cat# [TT210002]).