

Product datasheet for **TP325131M**

C1orf144 (SZRD1) (NM_001114600) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens chromosome 1 open reading frame 144 (C1orf144), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225131 representing NM_001114600 Red =Cloning site Green =Tags(s) MEDEEVAESWEEAADSGEIDRRLEKKLKITQKESRKS KSPPKVPVIVIQDDSLPAGPPPQIRILKRPTSNG VVSSPNSTSRPTLPVKSLAQREAEYAEARKRILGSASPEEEQEKPI LDRPTRISQPEDSRQPNNVIRQPL GPDGSQGFKQRR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	16.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001108072</u>
Locus ID:	26099
UniProt ID:	<u>Q7Z422</u>




[View online »](#)

Cytogenetics: 1p36.13

RefSeq ORF: 456

Synonyms: C1orf144

Product images:



Coomassie blue staining of purified SZRD1 protein (Cat# [TP325131]). The protein was produced from HEK293T cells transfected with SZRD1 cDNA clone (Cat# [RC225131]) using MegaTran 2.0 (Cat# [TT210002]).