

## Product datasheet for TP325198

### CBX1 (NM\_001127228) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromobox homolog 1 (HP1 beta homolog Drosophila ) (CBX1), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225198 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MGKKQNKKKVEEVLEEEEEYVVEKVLDRRWVKGKVEYLLKWKGFSDDEDNTWEPEENLDCPDLIAEFLQS QKTAHETDKSEGGKRKADSDSEDKGEESKPKKKKEESEKPRGFARGLEPERIIGATDSSGELMFLMKWKN SDEADLVPAKEANVKCPQVVISFYEERLTWHSYPSEDDDKKDDKN
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	21.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001120700</a></u>
Locus ID:	10951
UniProt ID:	<u><a href="#">P83916</a></u> , <u><a href="#">Q6IBN6</a></u>




[View online »](#)

RefSeq Size:	2253
Cytogenetics:	17q21.32
RefSeq ORF:	555
Synonyms:	CBX; HP1-BETA; HP1Hs-beta; HP1Hsbeta; M31; MOD1; p25beta

**Summary:** This gene encodes a highly conserved nonhistone protein, which is a member of the heterochromatin protein family. The protein is enriched in the heterochromatin and associated with centromeres. The protein has a single N-terminal chromodomain which can bind to histone proteins via methylated lysine residues, and a C-terminal chromo shadow-domain (CSD) which is responsible for the homodimerization and interaction with a number of chromatin-associated nonhistone proteins. The protein may play an important role in the epigenetic control of chromatin structure and gene expression. Several related pseudogenes are located on chromosomes 1, 3, and X. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified CBX1 protein (Cat# TP325198). The protein was produced from HEK293T cells transfected with CBX1 cDNA clone (Cat# [RC225198]) using MegaTran 2.0 (Cat# [TT210002]).