

Product datasheet for TP325206M

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PARK7 (NM_001123377) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human Parkinson disease (autosomal recessive, early onset) 7

(PARK7), transcript variant 2, 100 μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC225206 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

 $MASKRALVILAKGAEEMETVIPVDVMRRAGIKVTVAGLAGKDPVQCSRDVVICPDASLEDAKKEGPYDVV\\VLPGGNLGAQNLSESAAVKEILKEQENRKGLIAAICAGPTALLAHEIGFGSKVTTHPLAKDKMMNGGHYT$

YSENRVEKDGLILTSRGPGTSFEFALAIVEALNGKEVAAQVKAPLVLKD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 19.7 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001116849

Locus ID: 11315

UniProt ID: Q99497, V9HWC2





RefSeq Size: 921

Cytogenetics: 1p36.23 RefSeq ORF: 567

Synonyms: DJ-1; DJ1; GATD2; HEL-S-67p

Summary: The product of this gene belongs to the peptidase C56 family of proteins. It acts as a positive


regulator of androgen receptor-dependent transcription. It may also function as a redox-sensitive chaperone, as a sensor for oxidative stress, and it apparently protects neurons against oxidative stress and cell death. Defects in this gene are the cause of autosomal recessive early-onset Parkinson disease 7. Two transcript variants encoding the same protein

have been identified for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Protease

Protein Pathways: Parkinson's disease

Product images:



Coomassie blue staining of purified PARK7 protein (Cat# [TP325206]). The protein was produced from HEK293T cells transfected with PARK7 cDNA clone (Cat# [RC225206]) using MegaTran 2.0 (Cat# [TT210002]).